Steerable Electrophysiology Catheter

Instructions for use

Legal Manufacturer
CathRx Ltd
Unit 8, 2-8 South Street,
Rydalmere NSW 2116
Australia
Tel: + 61 2 7804 1846
Email: info@cathrx.com
Website: www.khelixmedical.com




European Authorised Representative
Advena Ltd
Tower Business Centre
2nd Flr, Tower Street
Swatar, BKR 4013
Malta

CONTENTS

Images	2
English	3
Warranty and Disclaimer	5
Deutsch (German)	6
Français (French)	8
Italiano (Italian)	10
Nederlands (Dutch)	12
Español (Spanish)	14
Português (Portuguese)	16
Svenska (Swedish)	18
Česky (Czech)	20
Dansk (Danish)	22
Norsk (Norwegian)	24
Symbol Definition	26

IMAGES

Steerable Electrophysiology Catheter

DEVICE DESCRIPTION

The Khelix Steerable Electrophysiology Catheter is designed to facilitate electrophysiological mapping of the heart. It consists of a handle and a shaft containing an array of platinum-iridium electrodes for stimulation and recording. The device is available in a number of different electrode configurations and curve types.

The distal tip can be deflected in one direction by pushing the deflection knob (A), as shown in [Image 1]. The tip is straightened by returning the knob to the original position.

INDICATIONS

The Khelix Steerable Electrophysiology Catheter is indicated for the electrical recording or stimulation of endocardial structures.

CONTRAINDICATIONS

- The Khelix Steerable Electrophysiology Catheter has not been shown to be safe and effective for radio frequency (RF) ablation.
- The transeptal approach is contraindicated in patients with left atrial thrombus, myxoma, or an inter-atrial patch or baffle.
- The use of this catheter may not be appropriate for use in patients with prosthetic valves.
- Electrophysiology procedures are contraindicated for patients with unstable cardiac conditions; e.g. acute myocardial infarction, unstable angina, hemodynamic instability.
- There is a relative contraindication for cardiac catheterisation procedures in patients with active systemic infection.

ADVERSE REACTIONS

A number of serious adverse reactions have been documented for cardiac catheterisation procedures, including cardiac tamponade, pulmonary embolism, myocardial infarction, stroke and death.

The following additional complications associated with cardiac catheterisation have also been reported in the literature: vascular bleeding, local haematoma, thrombosis, AV fistula, pseudoaneurysm, thromboembolism, vasovagal reaction, cardiac perforation, air embolism, arrhythmia, valvular damage, pneumothorax and haemothorax.

DIRECTIONS FOR USE

The Khelix Steerable Electrophysiology Catheter interfaces with standard electrophysiological recording equipment, such as ECG monitoring equipment, impedance based navigational equipment and cardiac stimulation equipment, via the CathRx Extension Cable (supplied separately).

Start of Procedure

1. Remove the catheter and external cable from their packaging using aseptic technique. Place the devices on a sterile working area. Ensure the catheter is not kinked.
2. Inspect the catheter and the extension cable for damage to the insulation or connectors. The connector contacts must be dry.
3. Test all functions of the catheter. When testing the deflection mechanism, ensure the catheter is held straight rather than coiled as shown in [Image 2].
4. Connect the extension cable by inserting the distal end of the cable into the socket located at the proximal end of the catheter handle. Prior to connection, ensure the key (raised ridges) on the cable connector (B) are aligned with the key slots inside the connector socket of the catheter handle. To prevent catheter shaft deflection during connection, grip the catheter deflection knob (C) during cable insertion [Image 3].
5. Ensure the catheter is un-deflected prior to insertion into the vasculature or into an introducer.

During Procedure

1. Advance the catheter through the vascular access site into the desired intracardiac position using fluoroscopic and ECG guidance. Vascular access is created with an appropriately sized haemostasis introducer sheath to accommodate the catheter.

At End of Procedure

1. Prior to removing the catheter, straighten the catheter by pulling the deflection knob to the neutral position.
2. Remove the catheter from the patient and dispose the catheter according to standard hospital procedure for biologically contaminated material.
3. Disconnect extension cable from catheter by gripping the extension cable pin head where the arrows are printed on the pin head and pull outwards (B). To prevent catheter shaft deflection during disconnection, grip the catheter deflection knob (C) during disconnection [Image 3].
4. Disconnect extension cable from system, immediately wipe down the cable with a sponge or

towel moistened with water. Refer extension cable Instructions for Use for cleaning and sterilisation instructions.

WARNINGS AND PRECAUTIONS

- Do not attempt to operate the catheter prior to completely reading and understanding these Instructions for Use.
- The catheter must only be used with the CathRx Extension Cable.
- The catheter should only be used by clinicians trained in cardiac electrophysiology procedures, in a fully equipped electrophysiology laboratory.
- Do not use if the package is open or damaged.
- Use the device prior to the 'Use By' date on the package label.
- The catheter is for single use only. Do not reuse, reprocess or re-sterilise the devices. Reuse, reprocessing or re-sterilisation may compromise the structural integrity of the devices and/or lead to device failure.
- Use of anti-coagulants should be considered for left-side procedures to reduce risk of thrombi formation and stroke.
- Cardiac catheterisation procedures present the potential for significant X-ray exposure. Exposure to X-rays can result in acute radiation injury as well as increased risk for somatic and genetic effects in both patients and electrophysiology laboratory staff, due to the X-ray beam intensity and duration of the fluoroscopic imaging. Steps should be taken to minimise X-ray exposure.
- In view of the potential for X-ray exposure and the associated risks, careful consideration must be given for the use of the device in pregnant women.
- Care should be taken when introducing, manipulating and withdrawing the catheter to avoid cardiac trauma and possible damage to the catheter which may result in detachment of catheter components.
- The catheter electrodes, the catheter and extension cable connectors, and the dispersive electrode, should not contact other conductive parts and earth.
- To avoid char formation on the catheter electrodes, do not apply RF energy when an ablation catheter is in contact with one or more of the Khelix Steerable Electrophysiology Catheter electrodes.
- Do not expose the catheter to organic solvents such as alcohol.
- Do not immerse the catheter handle or any cable connector in fluid as electrical performance could be affected.
- Do not use the catheter system near MRI equipment, as movement or heating of the catheter may occur and the image or display may become distorted.
- The catheter is MR unsafe.
- The catheter is a type CF applied part.
- Any equipment used in connection with the Khelix Steerable Electrophysiology Catheter and its accessories must be defibrillation-proof type CF, and meet IEC 60601-1 and IEC 60601-1-2.
- The catheter, together with its accessories has been tested to allow connection with stimulation equipment with a rated current of 25mA, and a rated voltage of 27V.
- Store in a cool, dry, dark place (minimum temperature: -10°C; maximum temperature: 55°C; maximum relative humidity: 90%).
- If the catheter becomes kinked, dispose of it. Replace with an undamaged device

WARRANTY AND DISCLAIMER

CathRx warrants that if it determines that a product was defective or faulty in materials and/or workmanship at the time of purchase, CathRx will, at its option, provide a replacement CathRx product or refund the purchase price of the defective or faulty product.

This limited warranty only applies if the following conditions are met:

- The product was packaged and labelled by CathRx;
- The product is returned to CathRx for evaluation within 30 days of identification of the defect or fault by the original purchaser;
- The product has not been repaired, altered, modified, mishandled or reprocessed inappropriately;
- The product has been used, stored, cleaned, sterilised and reprocessed in accordance with the product labelling and these Instructions for Use; and
- The product is not used after the 'Use By' date marked on the packaging.

OTHER THAN AS EXPRESSLY SET OUT IN THIS LIMITED WARRANTY, TO THE MAXIMUM EXTENT PERMITTED BY LAW, CATHRX MAKES NO WARRANTY IN RELATION TO THE PRODUCT, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION ANY WARRANTY OF SAFETY, MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR THAT THE PRODUCT WILL BE DEFECT OR FAULT FREE. TO THE MAXIMUM EXTENT PERMITTED BY LAW, CATHRX EXCLUDES ALL SUCH WARRANTIES.

CathRx accepts no liability and the purchaser of the product assumes all liability, whether based on warranty, tort, contract, negligence, under statute or otherwise, for any kind of loss or damage (including special, incidental or consequential) or injury (including death) arising directly or indirectly as a result of or in connection with the products, including arising as a result of the handling, possession, use or misuse of the product. CathRx's entire liability, including for a breach of warranty or condition which cannot lawfully be modified or excluded, is limited, at the option of CathRx, to either the replacement of the CathRx product or the refund of the purchase price of the product. The remedies set forth in this limited warranty are the exclusive remedy available to any person. The use of the product is deemed to be acceptance of the terms and conditions of this limited warranty.

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