## Diagnostic Extension Cable

|                         |         |        |
|-------------------------|---------|--------|
| Instructions for Use    | English | Page 2 |
| Warranty and Disclaimer | English | Page 7 |

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European Authorised  
Representative  
Advena Ltd  
Tower Business Centre  
2nd Flr, Tower Street  
Swatar, BKR 4013  
Malta

## DESCRIPTION AND INTENDED PURPOSE

The Khelix Diagnostic Extension Cable provides the interface between a Khelix diagnostic electrophysiology catheters and standard electrophysiological recording equipment, such as ECG monitoring equipment, impedance based navigational equipment and cardiac stimulation equipment. Khelix diagnostic extension cable are available in both 10 pin for EP catheters with 4 to 10 electrodes and 20 pin for EP catheters with more than 10 electrodes. An example of 20 pin diagnostic extension cable is shown below.


Image: Khelix diagnostic extension cable 20 pin.

The cable is provided by the manufacturer sterile. It may re-used subject to the cleaning and sterilisation instructions below.

## CONTRAINDICATIONS


There are no known contraindications to this cable.

## DIRECTIONS FOR USE

### Start of Procedure


- Step 1 Connect the extension cable by inserting the distal end of the cable into the socket located at the proximal end of the catheter handle. Prior to connection, ensure the key (raised ridges) on the cable connector are aligned with the key slots inside the

connector socket of the catheter handle. To prevent catheter shaft deflection during connection, grip the catheter deflection knob during cable insertion.

### End of Procedure

**Step 1** Disconnect extension cable from catheter by gripping the extension cable pin head where the arrows are printed on the pin head and pull outward. To prevent catheter shaft deflection during disconnection, grip the catheter deflection knob during disconnection.

Step 2 Disconnect extension cable from system, immediately wipe down the cable with a sponge or towel moistened with water. Refer extension cable Instructions for Use for cleaning and sterilisation instructions.

## **CLEANING AND STERILISATION**

The extension cable may be sterilised using steam (autoclave) for reuse. Prior to autoclaving, the cable must first be cleaned.

The manufacturer has validated the cable for twenty (20) cleaning/sterilisation cycles, using the procedures described below.

### **Cleaning Procedure**

1. Completely disconnect the extension cable from the system.
2. Immediately wipe down the cable with a sponge or towel moistened with water.
3. Contain all contaminated cables to minimise the risk of exposing employees to pathogens and prevent damage to the device.
4. Cover the cable with a towel that has been moistened with water to prevent body fluids from drying.

#### **5. Manual Cleaning Procedure**

1. Rinse the cable in running tap water for 30 seconds.
2. Prepare a solution by diluting 30 mL of Enzol Enzymatic Detergent (or equivalent) in 3785mL of tap water, or by following the manufacturers' instructions.
3. Soak the extension cable in the enzymatic detergent solution for a minimum of 5 minutes.
4. Remove device from detergent solution and brush all surfaces for a minimum of 2 minutes using a non-abrasive nylon bristle (9.5mm) general instrument brush (Key Surgical Item N- 3000 or equivalent) saturated with the enzymic detergent solution.
5. Rinse the cable in running tap water for a minimum of 30 seconds.
6. Immerse the cable in deionised water 3 separate times for 10 seconds each.
7. Wipe the device with a non-linting wipe and allow to completely dry.

#### **6. Automatic Cleaning Procedure**

Use a Belimed Washer Model 250 (or equivalent) with an enzymatic detergent (Ecolab or equivalent). Run the cycle with the following paramters:

1. Wash 45°C 4 minutes dose pump 4 (detergent) 5 mL
2. Wash 60°C 3 minutes
3. Rinse with unheated water for 1 minute
4. Rinse 60°C 1 minute
5. Thermal Disinfection A<sub>0</sub> 93°C
6. A<sub>0</sub> value: A<sub>0</sub>3000
7. Dry 123°C air 14 minutes

**WARNING**

Failure to properly clean the device may lead to inadequate sterilisation.

**Sterilisation Procedure**

Place cable in a CathRx Khelix Sterilization Tray – P/N SPC0057 and wrap tray with an appropriate grade sterilization wrap designated for use with a tray (Kimgard One Step KC400 or equivalent).

Only use the following steam sterilisation autoclave parameters to sterilise the cable.

1. Cycle Type: Dynamic Air Removal
2. Pre-vacuum pressure: 0.914 bar
3. Pre-vacuum Pulses: 3
4. Dwell Temperature: 132°C
5. Dwell Pressure: 1.85 bar
6. Dwell Time: 4 minutes
7. Drying Time: 20 minutes

The instructions provided above have been validated by the manufacturer of the cable as being CAPABLE of preparing a medical device for re-use. It remains the responsibility of the processor to ensure that the reprocessing as actually performed using equipment, materials and personnel in the reprocessing facility achieve the desired result. This normally requires validation and routine monitoring of the process.


**Ethylene Oxide Sterilisation**

Testing by the manufacturer has demonstrated that the cable can withstand 20 cycles of ethylene oxide sterilisation. However, as ethylene oxide sterilisation methods require sterilisation equipment-specific validation of parameters, it is the responsibility of the user to validate the effectiveness of the procedure with their sterilisation equipment. This validation should include a method for determining sterility such as biological indicators.

## WARNINGS AND PRECAUTIONS

- Do not use if the package is open or damaged.
- Use the device prior to the 'Use By' date on the package label.
- Do not re-use the cable beyond the validated 20 cycles of cleaning/sterilisation.
- The extension cable connectors should not contact other conductive parts and earth.
- Any equipment used with the cable must be defibrillation-proof type CF, and meet IEC 60601-1 and IEC 60601-1-2.
- The cable, together with catheters, has been tested to allow connection with stimulation equipment with a rated current of 25mA, and a rated voltage of 27V.
- Store in a cool, dry, dark place (minimum temperature: -10°C; maximum temperature: 55°C; maximum relative humidity: 90%).
- If the cable becomes kinked or otherwise damaged, dispose of it. Replace with an undamaged device.

## SYMBOL DEFINITIONS

|                                                                                     |                                                     |                                                                                     |                                            |
|-------------------------------------------------------------------------------------|-----------------------------------------------------|-------------------------------------------------------------------------------------|--------------------------------------------|
|     | Item number for ordering                            |    | Humidity limitation                        |
|     | Sterilised with ethylene oxide gas                  |    | Keep dry                                   |
|     | Manufacturer                                        |    | Keep away from sunlight                    |
|     | Lot number                                          |    | MR unsafe                                  |
|   | Device reference number                             |  | Serial number                              |
|   | Authorised representative in the European Community |  | Refer to accompanying Instructions for Use |
|  | Use by                                              |  | Temperature limitation                     |
|   | Date of manufacture                                 |  | Do not use if package is damaged           |
|   | Defibrillator-proof Type CF Applied Part            |                                                                                     |                                            |

## **WARRANTY AND DISCLAIMER**

CathRx warrants that if it determines that a product was defective or faulty in materials and/or workmanship at the time of purchase, CathRx will, at its option, provide a replacement CathRx product or refund the purchase price of the defective or faulty product.

This limited warranty only applies if the following conditions are met:

- The product was packaged and labelled by CathRx;
- The product is returned to CathRx for evaluation within 30 days of identification of the defect or fault by the original purchaser;
- The product has not been repaired, altered, modified, mishandled or reprocessed inappropriately;
- The product has been used, stored, cleaned, sterilised and reprocessed in accordance with the product labelling and these Instructions for Use; and
- The product is not used after the 'Use By' date marked on the packaging.

OTHER THAN AS EXPRESSLY SET OUT IN THIS LIMITED WARRANTY, TO THE MAXIMUM EXTENT PERMITTED BY LAW, CATHRX MAKES NO WARRANTY IN RELATION TO THE PRODUCT, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION ANY WARRANTY OF SAFETY, MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR THAT THE PRODUCT WILL BE DEFECT OR FAULT FREE. TO THE MAXIMUM EXTENT PERMITTED BY LAW, CATHRX EXCLUDES ALL SUCH WARRANTIES.

CathRx accepts no liability and the purchaser of the product assumes all liability, whether based on warranty, tort, contract, negligence, under statute or otherwise, for any kind of loss or damage (including special, incidental or consequential) or injury (including death) arising directly or indirectly as a result of or in connection with the products, including arising as a result of the handling, possession, use or misuse of the product. CathRx's entire liability, including for a breach of warranty or condition which cannot lawfully be modified or excluded, is limited, at the option of CathRx, to either the replacement of the CathRx product or the refund of the purchase price of the product. The remedies set forth in this limited warranty are the exclusive remedy available to any person. The use of the product is deemed to be acceptance of the terms and conditions of this limited warranty.

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