# Steerable Electrophysiology Catheter with Reuse Functionality

Lenkbare Elektrophysiologie-Katheter Mit Wiederverwendungsfunktion | Cathéter Électrophysiologique Orientable Avec Fonctionnalité de Réutilisation | Catetere Orientabile Per Elettrofisiologia Con Funzionalità Di Riutilizzo | Stuurbare Elektrofysiologiekatheter Met Functionaliteit Voor Hergebruik | Catéter de Electrofisiologia Orientable Con Funcionalidad de Reutilización | Cateter de Eletrofisiologia Direcionável Com Funcionalidade de Reutilização | Styrbara Elektrofysiologikateter Med Återanvändningsfunktionalitet | Řiditelný Elektrofyzilogický Katetr S Vícenásobně Využitelnými Komponentami | Styrbare Elektrofysiologiske Kateter Med Genanvendelsesfunktionalitet | Styrbart Elektrofysiologikateter Med Gjenbruksfunksjonalitet

## Instructions for use

Gebruiksaanweisung | Conseils d'utilisation | Istruzioni d'uso | Gebruiksaanwijzingen | Instrucciones de uso | Instruções de uso | Bruksanvisning | Návod k použití | Brugsanvisning | Bruksanvisning

CathRx Ltd  
Unit 8/2-8 South Street  
Rydalmeere NSW 2116  
Australia  
Tel: (+61) 2 7804 1846  
Email: [info@cathrx.com](mailto:info@cathrx.com)  
[www.cathrx.com](http://www.cathrx.com)


Advena Ltd  
Tower Business Centre  
2<sup>nd</sup> Flr.  
Tower Street  
Swatar, BKR 4013  
Malta








[www.cathrx.com](http://www.cathrx.com)

## CONTENTS

Images	2
English	4
Deutsch (German)	9
Français (French)	15
Italiano (Italian)	21
Nederlands (Dutch)	27
Español (Spanish)	33
Português (Portuguese)	39
Svenska (Swedish)	45
Čeština (Czech)	51
Dansk (Danish)	56
Norsk (Norwegian)	61
Symbol Definition	66

## KHELIX STEERABLE ELECTROPHYSIOLOGY CATHETER WITH REUSE FUNCTIONALITY

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- STERILE. Sterilized with Ethylene Oxide.
- Do not use if the package is open or damaged.
- Handle/Guidewire is reusable per reprocessing instructions provided below.
- Catheter sheath is for single use only. Do not reuse, reprocess or re-sterilise.
- Use the device prior to the "Use By" date on the package label.

These instructions apply to the following products:

Table 1: Product Names	
Product Name	Product Reference Number (REF)
Steerable Electrophysiology Catheters with Reuse Functionality	MS1-D0U-04A6252
	MS1-D0U-10A6252
	MS1-F0U-10A6252
	MS2-P0U-20P62122
	MS1-G0U-10A6282
	MS1-F0U-10A6282
	MS1-G0U-10A6252
	MS2-Q0U-20P62122
Reusable Handle/Guidewire*	MS1-D0U
	MS1-F0U
	MS2-P0U
Catheter Sheath Single Use	MS1-04A6252
	MS1-10A6252
	MS2-20P62122
	MS1-10A6282

\* Reusable Handle/Guidewire is referred similarly as Steerable Electrophysiology Catheter with Reuse Functionality or reusable handle/guidewire in this IFU.

**Sterilization:** This product and its packaging have been sterilized with ethylene oxide gas (EO). Even though the product is processed in compliance with all applicable laws and regulations relating to EO exposure, Proposition 65, a State of California voter initiative, requires the following notice:

Warning: This product and its packaging have been sterilized with ethylene oxide. The packaging may expose you to ethylene oxide, a chemical known to the State of California to cause cancer or birth defects or other reproductive harm.

### DEVICE DESCRIPTION

Khelix Steerable Electrophysiology Catheter with Reuse Functionality is designed to facilitate electrophysiological mapping of the heart. It consists of a handle (A) and catheter shaft (B) [Image 1].

The catheter shaft contains an array of platinum-iridium electrodes for stimulation and recording, and an internal guidewire.

The handle has a deflection knob (C) [Image 1]. Push the deflection knob away from the handle to deflect the distal tip of the catheter sheath via the guidewire. To straighten the tip, return the knob to the neutral position. The device is provided sterile and fully assembled for use. After use it can be disassembled into:

1. Catheter sheath (D) [Image 2], and
2. Catheter handle (E) / guidewire (F) [Image 2].

The handle/guidewire may be reused after cleaning and sterilisation. The sheath is for single use only and must be discarded after use. Both components are also available and sold separately.

### Catheter Sheath Single Use

Available in a number of different electrode configurations. Intended to be assembled with a sterile handle/guidewire to form a Steerable Electrophysiology Catheter with Reuse Functionality.

### Reusable Handle/Guidewire

Available with either a 10-pin or a 20-pin electrical connector (for use with catheter sheaths as per Table 2), and in different guidewire curve types. Intended to be assembled with a sterile catheter sheath to form a Steerable Electrophysiology Catheter with Reuse Functionality.

Table 2: Pin Connector Types		
Handle/ Guidewire	REF prefix	Catheter sheath
10-pin connector	MS1	4–10 electrodes
20-pin connector	MS2	12–20 electrodes

Khelix Steerable Electrophysiology Catheter with Reuse Functionality interfaces with standard electrophysiological recording equipment, such as ECG monitoring equipment, impedance based navigational equipment and cardiac stimulation equipment, via the Khelix Diagnostic Extension Cable (supplied separately).

Khelix Steerable Electrophysiology Catheter with Reuse Functionality is intended to be used with the following accessories:

Table 3: Accessories for Khelix Steerable Electrophysiology Catheter with Reuse Functionality
Khelix Diagnostic Extension Cable
Khelix Sterilisation Tray (Ref: A-T01)

Both accessories are available and sold separately by CathRx Ltd

## INDICATIONS

Khelix Steerable Electrophysiology Catheter with Reuse Functionality is indicated for the electrical recording or stimulation of endocardial structures. The handle/guidewire may be reused after cleaning, sterilisation and re-assembly with a new catheter sheath.

## INTENDED USER

Khelix Steerable Electrophysiology Catheter with Reuse Functionality is used by interventional cardiologists in electrophysiology (EP) studies to determine the cause of an abnormal heart rhythm, to locate the site of origin of an abnormal heart rhythm, to decide the best treatment for an abnormal heart rhythm, and/or to check the effectiveness of ablation therapy.

## CONTRAINDICATIONS

- The Khelix Steerable Electrophysiology Catheter with Reuse Functionality has not been shown to be safe and effective for radio frequency (RF) ablation.
- The transeptal approach is contraindicated in patients with left atrial thrombus, myxoma, or an inter-atrial patch or baffle.
- The use of this catheter may not be appropriate for use in patients with prosthetic valves.
- Electrophysiology procedures are contraindicated for patients with unstable cardiac conditions, e.g. acute myocardial infarction, unstable angina, hemodynamic instability.
- There is a relative contraindication for cardiac catheterisation procedures in patients with active systemic infection.

## ADVERSE REACTIONS

A number of serious adverse reactions have been documented for cardiac catheterisation procedures, including cardiac tamponade, pulmonary embolism, myocardial infarction, stroke and death.

The following additional complications associated with cardiac catheterisation have also been reported in the literature: vascular bleeding, local haematoma, thrombosis, AV fistula, pseudoaneurysm, thromboembolism, vasovagal reaction, cardiac perforation, air embolism, arrhythmia, valvular damage, pneumothorax and haemothorax.

## DIRECTIONS FOR USE

### At point of use

1. Inspect the catheter and packaging before opening. The contents of the package are sterile unless the package is opened or damaged.

2. If the package is opened or damaged, or if the package was opened and the catheter unused, do not use the catheter.
3. A detachable package label has been provided and may be affixed to the medical record of the patient.

### Start of Procedure

1. Remove new catheter from its sterile packaging using aseptic technique or assemble the sterilised or new reusable handle/guidewire and new catheter sheath according to assembly instructions provided below.
2. Place the catheter on a sterile working area.
3. Inspect the catheter for any damages. Do not use the catheter if kinked.
4. Inspect the catheter for damage to the insulation or connectors. The connector contacts must be dry.
5. Test all functions of the catheter. When testing the deflection mechanism, ensure the catheter is straight and not coiled **[Image 3]**.
6. Connect the extension cable. Grip the deflection knob to prevent catheter shaft deflection **[Image 4]**. Align the key (raised ridges) on the cable connector **(G)** with the key slots inside the socket of the reusable handle/guidewire and insert the connector.
 

 Khelix Steerable Electrophysiology Catheter with Reuse Functionality must only be used with Khelix Diagnostic Extension Cable.
7. Ensure the catheter is un-deflected prior to insertion into the vasculature or into an introducer.

### During Procedure

1. Vascular access is created with a minimum 6 French haemostasis introducer sheath to accommodate the catheter.
2. Advance the catheter through the vascular access site into the desired intracardiac position using fluoroscopic and ECG guidance.

### End of Procedure

1. Straighten the catheter by pulling the deflection knob to the neutral position.
2. Remove the catheter from the patient.
3. Grip the deflection knob to prevent catheter shaft deflection. Disconnect extension cable from handle, gripping the extension cable pin head at the arrows and pulling outwards **[Image 5]**.
4. If the handle/guidewire will not be reused, dispose of the complete catheter assembly according to standard hospital procedure for biologically contaminated material.

## DISASSEMBY INSTRUCTIONS

1. Ensure the catheter is fully un-deflected, straight and disconnected from connecting cable before disassembly **[Image 3]**.

- Grip the deflection knob to prevent catheter shaft deflection. Gripping the deflection knob will ensure the guidewire remains un-deflected during sheath removal.

Deflection of the guidewire during sheath removal may result in difficulty in removing the sheath and potential damage to the guidewire.

- With the other hand pinch the blue collar **[Image 6]** and pull the sheath away from the handle until the sheath electrical connector **(H)** is clear **[Image 7]**.
- Smoothly retract the guidewire from the sheath, ensuring the sheath remains straight throughout. As the sheath is retracted, the guidewire may be coiled in one large loop.
- Once it is completely removed from the guidewire, dispose of the sheath according to standard hospital procedure for biologically contaminated material.
- Push deflection knob forward and wipe down all available surfaces of the reusable handle/guidewire with a lint-free cloth dampened with tap water. Pull deflection knob back to neutral position.
- Coil the guidewire in a single large loop and place the handle/guidewire in Khelix Sterilisation Tray A-T01 (or similar) **[Image 8]**.
- Clean, sterilise and transport the reusable handle/guidewire according to instructions provided below.

## REPROCESSING INSTRUCTIONS

Only the handle/guidewire assembly can be reprocessed following the instruction below, provided that the limitations for reprocessing indicated below is not exceeded.

### Cleaning (Manual)

- Prepare ENZOL® Enzymatic Detergent solution (or equivalent) in a container by diluting 30 ml of detergent in 3785 ml of tap water (20–25°C).

Wear gloves during solution preparation and cleaning procedures.

- Rinse device under running tap water (20–25°C) for 30 seconds.
- Submerge device in the cleaning solution for 5 minutes.
- Remove device from the cleaning solution and place on a clean CSR (sterilisation) drape.
- Using a non-abrasive nylon bristled brush saturated with the ENZOL® Enzymatic Detergent solution (or equivalent), gently brush all surfaces of the device for 2 minutes, including the surface area when the deflection knob has been pushed forward.
- Rinse device in running tap water (20–25°C) for 30 seconds.
- Immerse device in de-ionized (DI) water 3 separate times for 10 seconds each.

- Wipe the device using a non-linting wipe, including the surface area when the deflection knob has been pushed forward. Place on a clean CSR (sterilisation) drape and allow to dry until visibly free of moisture.
- After drying, inspect visually for evidence of soil remaining on device.

### Cleaning (Automated)

- Push deflection knob forward and wipe device using a non-linting wipe. Pull deflection knob back to neutral position and process the handle/guidewire assembly in a Belimed Washer Model WD250 using Ecolab® Neutral Enzymatic Detergent (or an equivalent automated washer/enzymatic detergent).
- Place handle/guidewire assembly in the washer equipment manufacturer's tray or a suitable CSD washing tray. Ensure that the guidewire does not kink.
- Washing parameters:
  - Wash at 45°C for 4 minutes. Dose pump 4 (detergent) 5 MI
  - Wash at 60°C for 3 minutes
  - Rinse with unheated water for 1 minute
  - Rinse at 60°C for 1 minute
  - Thermal Disinfection A0 93°C 2.5 minutes
  - A0 value: A03000
  - Dry at 123°C air for 14 minutes.
- Upon completion of wash cycle, use gloves to remove device(s). Wipe using a non-linting wipe. Place on a clean CSR (sterilisation) drape. Dry until visibly free of moisture.
- After drying, inspect visually for evidence of soil remaining on each device.

Failure to properly clean the device could lead to inadequate sterilisation.

### Inspection and Maintenance

- Visually inspect the reusable handle/guidewire for any damages.
- Test all functions of the catheter. When testing the deflection mechanism, ensure the catheter is straight and not coiled **[Image 3]**.
- If reusable handle/guidewire is not damaged or kinked, continue the use of reusable handle/guidewire and proceed to sterilisation.
- Dispose any damaged reusable handle/guidewire according to disposal instructions provided below.

### Sterilisation

- Place handle/guidewire assembly in a Khelix Sterilisation Tray A-T01 (Sold separately at CathRx Ltd) **[Image 8]**. Secure handle and guidewire in holders as indicated by the markings within the tray, ensuring the guidewire does not kink. Handle/guidewire assemblies are limited to one per tray.

2. Wrap the trays in a sterilisation wrap (eg Kimberly-Clark Kimguard One-Step Sterilization Wraps) using a simultaneous double-wrapping envelope fold technique.
3. The following steam sterilisation parameters were validated in an AMSCO Vacumatic "A" Autoclave:
  - Cycle Type: Dynamic Air Removal
  - Pre-vacuum pressure: 27" Hg (0.914 bar)
  - Pre-vacuum Pulses: 3
  - Dwell Temperature: 132°C
  - Dwell Pressure: 26.8 psig (1.85 bar)
  - Dwell Time: 4 minutes
  - Drying Time: 20 minutes.

## ASSEMBLY INSTRUCTIONS

The complete assembly procedure must be performed aseptically. No part of the sterile sheath or the sterile handle/guidewire must contact any object outside the sterile field.

1. Match the reusable handle/guidewire to the required catheter sheath (see **Table 2**). The deflection knob [**Image 9**] indicates the curve type (I) and the pin connector (J) of the reusable handle/guidewire.
2. Prepare an adequate sterile flat surface (e.g. sterile draped table) in the sterile field for the sheath and reusable handle/guidewire.
3. Assemble immediately prior to use.
4. Remove the reusable handle/guidewire from its sterile packaging or the sterilised wrapped tray using standard aseptic technique.
5. Check that the guidewire is fully un-deflected and straight, and ensure it is maintained in this position throughout the assembly process.

Always ensure the deflection knob is in the neutral position (un-deflected) during assembly and disassembly to avoid damage to the guidewire.

6. Place the handle/guidewire on the flat sterile surface, with the guidewire coiled in one large loop [**Image 10**].
7. Remove a new sheath from its sterile packaging and lay it straight on the flat sterile surface.
8. Hold handle and distal end of guidewire in one hand with the index finger and thumb approximately 2 cm from the guidewire tip [**Image 11**].
9. Hold the sheath electrical connector with the other hand. Insert guidewire into centre hole of connector and gently uncoil the guidewire to smoothly slide it all the way into the sheath [**Image 12**].
10. Insert the sheath electrical connector into the distal end of the handle and gently rotate until it is fully inserted and has clicked firmly into place [**Image 13**].
11. Check the catheter shaft can be smoothly deflected and straightened. Ensure that the catheter shaft is not coiled when this check is performed.

## LIMITATIONS ON REPROCESSING

Reprocessing Khelix Steerable Electrophysiology Catheter with Reuse Functionality has been validated for 10 cycles. One cycle consists of cleaning and sterilisation processes.

## STORAGE

Store Khelix Steerable Electrophysiology Catheter with Reuse Functionality in a cool, dry place (minimum temperature: -10°C [14°F]; maximum temperature: 55°C [131°F]; relative humidity: 85% ± 5%).

## CONTAINMENT AND TRANSPORT

1. Remove the used catheter sheath and dispose. Refer to disassembly instructions provided above.
2. Wipe the used reusable handle/guidewire with a lint-free cloth dampened with tap water to contain contaminants and prevent body fluids from drying.
3. Place and secure the used reusable handle/guidewire inside Khelix Sterilisation Tray (or similar). Place tray lid on and lock the tray.
4. Place the tray into an enclosed transport bin or enclosed case cart to prevent cross-contamination during transport.

## DISPOSAL

Dispose Khelix Steerable Electrophysiology Catheter with Reuse Functionality according to standard hospital procedure for biologically contaminated material.

## WARNINGS AND PRECAUTIONS

- Do not attempt to operate the catheter prior to completely reading and understanding these Instructions for Use.
- Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement and placement should be done under fluoroscopic guidance.
- Do not use excessive force to advance or withdraw the catheter when resistance is encountered.
- The catheter should only be used by clinicians trained in cardiac electrophysiology procedures, in a fully equipped electrophysiology laboratory.
- The catheter sheath is for single use only. Do not reuse, reprocess or resterilise the sheath. Reuse, reprocessing or resterilisation may compromise the structural integrity of the device and/or lead to device failure. For example, detachment of small components like tip, ring electrode dislocates, unable to sense signals in heart. Risk of cross contamination, poor catheter handling.
- Use of anti-coagulants should be considered for left-side procedures to reduce risk of thrombi formation and stroke.

- Cardiac catheterisation procedures present the potential for significant X-ray exposure. Exposure to X-rays can result in acute radiation injury as well as increased risk for somatic and genetic effects in both patients and electrophysiology laboratory staff, due to the X-ray beam intensity and duration of the fluoroscopic imaging. Steps should be taken to minimise X-ray exposure.
- In view of the potential for X-ray exposure and the associated risks, careful consideration must be given for the use of the device in pregnant women.
- Care should be taken when introducing, manipulating and withdrawing the catheter to avoid cardiac trauma and possible damage to the catheter which may result in detachment of catheter components.
- The catheter electrodes, the catheter and extension cable connectors, and the dispersive electrode, should not contact other conductive parts and earth.
- To avoid char formation on the catheter electrodes, do not apply RF energy when an ablation catheter is in contact with one or more of the Khelix Steerable Electrophysiology Catheter with Reuse Functionality electrodes.
- Do not expose the catheter to organic solvents such as alcohol.
- When connected to power source, do not immerse the catheter handle or any cable connector in fluid as electrical performance could be affected.
- The catheter is MR unsafe. Do not use the catheter system near MRI equipment, as movement or heating of the catheter may occur, and the image or display may become distorted.
- The catheter is a type CF applied part. Any equipment used in connection with this device and its accessories must be defibrillation-proof type CF and meet IEC 60601-1 and IEC 60601-1-2.
- The catheter, together with its accessories, has been tested to allow connection with stimulation equipment with a rated current of 25mA, and a rated voltage of 27V.
- If the catheter becomes kinked, dispose of it. Replace with an undamaged device.
- Serious incidents related to the device should be reported to the manufacturer and the competent authority of your country.

## WARRANTY AND DISCLAIMER

CathRx warrants that if it determines that a product was defective or faulty in materials and/or workmanship at the time of purchase, CathRx will, at its option, provide a replacement CathRx product or refund the purchase price of the defective or faulty product.

This limited warranty only applies if the following conditions are met:

- The product was packaged and labelled by CathRx;
- The product is returned to CathRx for evaluation within 30 days of identification of the defect or fault by the original purchaser;
- The product has not been repaired, altered, modified, mishandled or reprocessed inappropriately;
- The product has been used, stored, cleaned, sterilised and reprocessed in accordance with the product labelling and these Instructions for Use; and
- The product is not used after the 'Use By' date marked on the packaging.

OTHER THAN AS EXPRESSLY SET OUT IN THIS LIMITED WARRANTY, TO THE MAXIMUM EXTENT PERMITTED BY LAW, CATHRX MAKES NO WARRANTY IN RELATION TO THE PRODUCT, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION ANY WARRANTY OF SAFETY, MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR THAT THE PRODUCT WILL BE DEFECT OR FAULT FREE. TO THE MAXIMUM EXTENT PERMITTED BY LAW, CATHRX EXCLUDES ALL SUCH WARRANTIES.

CathRx accepts no liability and the purchaser of the product assumes all liability, whether based on warranty, tort, contract, negligence, under statute or otherwise, for any kind of loss or damage (including special, incidental or consequential) or injury (including death) arising directly or indirectly as a result of or in connection with the products, including arising as a result of the handling, possession, use or misuse of the product. CathRx's entire liability, including for a breach of warranty or condition which cannot lawfully be modified or excluded, is limited, at the option of CathRx, to either the replacement of the CathRx product or the refund of the purchase price of the product. The remedies set forth in this limited warranty are the exclusive remedy available to any person. The use of the product is deemed to be acceptance of the terms and conditions of this limited warranty.

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