Steerable Electrophysiology Catheter with Reuse Functionality

Lenkbare Elektrophysiologie-Katheter Mit Wiederverwendungsfunktion | Cathéter Électrophysiologique Orientable Avec Fonctionnalité de Réutilisation | Catetere Orientabile Per Elettrofisiologia Con Funzionalità Di Riutilizzo | Stuurbare Elektrofysiologiekatheter Met Functionaliteit Voor Hergebruik | Catéter de Electrofisiología Orientable Con Funcionalidad de Reutilización | Cateter de Eletrofisiologia Direcionável Com Funcionalidade de Reutilização | Styrbara Elektrofysiologikateter Med Återanvändningsfunktionalitet | Říditelný Elektrofyzilogický Katetr S Vícenásobně Využitelnými Komponentami | Styrbare Elektrofysiologiske Kateter Med Genanvendelsesfunktionalitet | Styrbart Elektrofysiologikateter Med Gjenbruksfunksjonalitet

Instructions for use

Gebruiksaanweisung | Conseils d'utilisation | Istruzioni d'uso | Gebruiksaanwijzingen | Instrucciones de uso | Instruções de uso | Bruksanvisning | Návod k použití | Brugsanvisning | Bruksanvisning

CathRx Ltd
Unit 8/2-8 South Street
Rydalmere NSW 2116
Australia
Tel: (+61) 2 7804 1846
Email: info@cathrx.com
www.cathrx.com


Advena Ltd
Tower Business Centre
2nd Flr.
Tower Street
Swatar, BKR 4013
Malta








www.cathrx.com

CONTENTS

Images	2
English	4
Deutsch (German)	10
Français (French)	17
Italiano (Italian)	24
Nederlands (Dutch)	31
Español (Spanish)	38
Português (Portuguese)	44
Svenska (Swedish)	51
Čeština (Czech)	57
Dansk (Danish)	63
Norsk (Norwegian)	69
Symbol Definition	75

KHELIX STEERABLE ELECTROPHYSIOLOGY CATHETER WITH REUSE FUNCTIONALITY

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- STERILE. Sterilized with Ethylene Oxide.
- Do not use if the package is open or damaged.
- Handle/Guidewire is reusable per reprocessing instructions provided below.
- The catheter sheath is for single use only. Do not reuse, reprocess or re-sterilise the sheath.
- Use the device prior to the "Use By" date on the package label.

These instructions apply to the following products:

Table 1: Product Names	
Product Name	Product Reference Number (REF)
Steerable Electrophysiology Catheters with Reuse Functionality	MS1-D0U-04A6252
	MS1-D0U-10A6252
	MS1-F0U-10A6252
	MS2-P0U-20P62122
	MS1-G0U-10A6282
	MS1-F0U-10A6282
	MS1-G0U-10A6252
	MS2-Q0U-20P62122
Reusable Handle/Guidewire*	MS1-D0U
	MS1-F0U
	MS2-P0U
Catheter Sheath Single Use	MS1-04A6252
	MS1-10A6252
	MS2-20P62122
	MS1-10A6282

* Reusable Handle/Guidewire is referred similarly as Steerable Electrophysiology Catheter with Reuse Functionality or handle/guidewire in this IFU.

Sterilization: This product and its packaging have been sterilized with ethylene oxide gas (EO). Even though the product is processed in compliance with all applicable laws and regulations relating to EO exposure, Proposition 65, a State of California voter initiative, requires the following notice:

Warning: This product and its packaging have been sterilized with ethylene oxide. The packaging may expose you to ethylene oxide, a chemical known to the State of California to cause cancer or birth defects or other reproductive harm.

DEVICE DESCRIPTION

Khelix Steerable Electrophysiology Catheter with Reuse Functionality is designed to facilitate electrophysiological mapping of the heart. It consists of a handle (A) and catheter shaft (B) [Image 1].

The catheter shaft contains an array of platinum-iridium electrodes for stimulation and recording, and an internal guidewire.

The handle has a deflection knob (C) [Image 1]. Push the deflection knob away from the handle to deflect the distal tip of the catheter sheath via the guidewire. To straighten the tip, return the knob to the neutral position. The device is provided sterile and fully assembled for use. After use it can be disassembled into:

1. Catheter sheath (D) [Image 2], and
2. Catheter handle (E) / guidewire (F) [Image 2].

The handle/guidewire may be reused after cleaning and sterilisation. The sheath is for single use only and must be discarded after use. Both components are also available separately.

Catheter Sheath Single Use

Available in a number of different electrode configurations. Intended to be assembled with a sterile handle/guidewire to form a Steerable Electrophysiology Catheter.

Reusable Handle/Guidewire

Available with either a 10-pin or a 20-pin electrical connector (for use with catheter sheaths per Table 2), and in different guidewire curve types. Intended to be assembled with a sterile catheter sheath to form a Steerable Electrophysiology Catheter.

Table 2: Pin Connector Types		
Handle/ Guidewire	REF prefix	Catheter sheath
10-pin connector	MS1	4–10 electrodes
20-pin connector	MS2	12–20 electrodes

The Khelix Steerable Electrophysiology Catheter with Reuse Functionality interfaces with standard electrophysiological recording equipment, such as ECG monitoring equipment, impedance based navigational equipment and cardiac stimulation equipment, via the Khelix Diagnostic Extension Cable (supplied separately).

Khelix Steerable Electrophysiology Catheter with Reuse Functionality is intended to be used with the following accessories:

Table 3: Accessories for Khelix Steerable Electrophysiology Catheter with Reuse Functionality
Khelix Diagnostic Extension Cable
Khelix Sterilisation Tray (Ref: A-T02)

Both accessories are available and sold separately by CathRx Ltd

INDICATIONS

Khelix Steerable Electrophysiology Catheter with Reuse Functionality is indicated for the electrical recording or stimulation of endocardial structures. The handle/guidewire may be reused after cleaning, sterilisation and re-assembly with a new catheter sheath.

INTENDED USER

Khelix Steerable Electrophysiology Catheter with Reuse Functionality is used by interventional cardiologists in electrophysiology (EP) studies to determine the cause of an abnormal heart rhythm, to locate the site of origin of an abnormal heart rhythm, to decide the best treatment for an abnormal heart rhythm, and/or to check the effectiveness of ablation therapy.

CONTRAINDICATIONS

- The Khelix Steerable Electrophysiology Catheter with Reuse Functionality has not been shown to be safe and effective for radio frequency (RF) ablation.
- The transeptal approach is contraindicated in patients with left atrial thrombus, myxoma, or an inter-atrial patch or baffle.
- The use of this catheter may not be appropriate for use in patients with prosthetic valves.
- Electrophysiology procedures are contraindicated for patients with unstable cardiac conditions, e.g. acute myocardial infarction, unstable angina, hemodynamic instability.
- There is a relative contraindication for cardiac catheterisation procedures in patients with active systemic infection.

ADVERSE REACTIONS

A number of serious adverse reactions have been documented for cardiac catheterisation procedures, including cardiac tamponade, pulmonary embolism, myocardial infarction, stroke and death.

The following additional complications associated with cardiac catheterisation have also been reported in the literature: vascular bleeding, local haematoma, thrombosis, AV fistula, pseudoaneurysm, thromboembolism, vasovagal reaction, cardiac perforation, air embolism, arrhythmia, valvular damage, pneumothorax and haemothorax.

DIRECTIONS FOR USE

At point of use

1. Inspect the catheter and packaging before opening. The contents of the package are sterile unless the package is opened or damaged.

2. If the package is opened or damaged, or if the package was opened and the catheter unused, do not use the catheter.
3. A detachable package label has been provided and may be affixed to the medical record of the patient.

Start of Procedure

1. Remove new catheter from its sterile packaging or assemble the sterilised or new reusable handle/guidewire and new catheter sheath according to assembly instructions provided below using aseptic technique. Place devices on a sterile working area.
2. Inspect the catheter for any damages. Do not use the catheter if kinked.
3. Inspect the catheter for damage to the insulation or connectors. The connector contacts must be dry.
4. Test all functions of the catheter. When testing the deflection mechanism, ensure the catheter is straight and not coiled **[Image 3]**.
5. Connect the extension cable. Grip the deflection knob to prevent catheter shaft deflection **[Image 4]**. Align the key (raised ridges) on the cable connector **(G)** with the key slots inside the socket of the handle and insert the connector.

The catheter must only be used with Khelix Diagnostic Extension Cable.

6. Ensure the catheter is un-deflected prior to insertion into the vasculature or into an introducer.

During Procedure

1. Vascular access is created with a minimum 6 French haemostasis introducer sheath to accommodate the catheter.
2. Advance the catheter through the vascular access site into the desired intracardiac position using fluoroscopic and ECG guidance.

End of Procedure

1. Straighten the catheter by pulling the deflection knob to the neutral position.
2. Remove the catheter from the patient.
3. Grip the deflection knob to prevent catheter shaft deflection. Disconnect extension cable from handle by gripping the extension cable pin head at the arrows and pulling outwards **[Image 5]**.
4. If not reusing the handle/guidewire assembly, dispose of the complete catheter assembly according to standard hospital procedure for biologically contaminated material.

DISASSEMBLY INSTRUCTIONS

1. Ensure the catheter is fully un-deflected, straight and disconnected from the extension cable before disassembly **[Image 3]**.

- Grip the deflection knob to prevent catheter shaft deflection. Gripping the deflection knob will ensure the guidewire remains un-deflected during sheath removal.
 -  Deflection of the guidewire during sheath removal may result in difficulty in removing the sheath and potential damage to the guidewire.
- With the other hand pinch the blue collar **[Image 6]** and pull the sheath away from the handle until the sheath electrical connector **(H)** is clear **[Image 7]**.
- Smoothly retract the guidewire from the sheath, ensuring the sheath remains straight throughout. As the sheath is retracted, the guidewire may be coiled in one large loop.
- Once it is completely removed from the guidewire, dispose of the sheath according to standard hospital procedure for biologically contaminated material.
- Push deflection knob forward and wipe down all available surfaces of handle/guidewire with a lint-free cloth dampened with tap water. Pull deflection knob back to neutral position.
- Coil the guidewire in a single large loop and place the handle/guidewire in Khexil Sterilisation Tray A-T02 (or similar) **[Image 8]**.
- Clean, sterilise and transport the reusable handle/guidewire according to instructions provided below.

REPROCESSING INSTRUCTIONS

-  Only the handle/guidewire assembly can be reprocessed following the instructions below, provided that the limitations for reprocessing indicated below is not exceeded.

Manual Cleaning Method

- Prepare ENZOL® Enzymatic Detergent solution (or equivalent) in a container by diluting 60 ml of ENZOL® Enzymatic Detergent in a gallon (3785 ml) of tap water at a temperature of 20–25°C.
 -  Wear gloves during solution preparation and cleaning procedures.
- Rinse device under running tap water (20–25°C) for 30 seconds.
- Soak device in the detergent solution for 5 minutes.
- Remove device from the detergent solution and place on a clean CSR (sterilisation) drape.
- Using a non-abrasive nylon bristled brush saturated with the detergent solution, gently brush all surfaces of the device for 2 minutes, including the surface area when the deflection knob has been pushed forward.
- Rinse device in running tap water (20–25°C) for 30 seconds.

- Immerse device in RO/DI water 3 separate times for 10 seconds each, using a new bath of RO/DI water each time.
- Wipe the device using a non-linting wipe, including the surface area when the deflection knob has been pushed forward.
- Place device on a clean CSR (sterilisation) drape and allow to dry until visibly free of moisture.
- After drying, inspect visually for evidence of soil remaining on device. If visible soil is present, remove soil by brushing using soft sponge/cloth or soft-bristled brush, and/or repeat of manual cleaning prior to sending to sterilisation.

Automated Cleaning Method

- Place the handle/guidewire in the Khexil Sterilisation Tray A-T02 (or similar) **[Image 8]**. Secure the handle/guidewire in the support brackets, ensuring the guidewire does not kink. Handle/guidewire assemblies are limited to one per tray.
- Check that the reusable handle/guidewire is properly fitted into the support brackets and is not damaging the device.
- Place tray lid on and lock the tray.
- Place the tray with the handle/guidewire in a medical grade washer (e.g. Steris AMSCO models or Belimed WD models).
- Start the automated cleaning using the validated washing parameters below:

	Detergent Type	Neutral enzymatic detergent
Wash 1	Water Quality	Tap Water
	Volume (dose pump)	¼ oz per gallon (1.87 g/l)
	Temperature	45°C tap water (set point)
	Minimum Time	4 minutes
Wash 2	Detergent Type	Neutral enzymatic detergent
	Water Quality	Tap Water
	Volume (dose pump)	¼ oz per gallon (1.87 g/l)
	Temperature	60°C tap water (set point)
Rinse	Minimum Time	3 minutes
	Water Quality	Tap Water
	Temperature	60°C tap water (set point)
Dry Time	Minimum Time	1 minute
	Temperature	115°C (set point)
	Minimum Time	14 minutes

- Upon completion of automated cleaning, while wearing gloves, remove the handle/guidewire from the tray.

7. Wipe the device using a non-linting wipe. Place on a clean CSR (sterilisation) drape and allow to dry until visibly free of moisture.
8. After drying, inspect the visually for evidence of soil remaining on the device. If visible soil is present, remove soil by manual cleaning and brushing using soft sponge/cloth or soft-bristled brush, and/or repeat of automatic cycles prior to sending to sterilisation.

Failure to properly clean the device could lead to inadequate sterilisation.

Inspection and Maintenance

1. Visually inspect the handle/guidewire for any damages.
2. Test all functions of the catheter. When testing the deflection mechanism, ensure the catheter is straight and not coiled **[Image 3]**.
3. If handle/guidewire is not damaged or kinked, continue the use of handle/guidewire and proceed to sterilisation.
4. Dispose any damaged handle/guidewire according to disposal instructions provided below.

Steam Sterilisation

Note: skip step 1 to 3 below if the reusable handle/guidewire is already placed into the tray during automated washing cycle.

1. Place the cleaned handle/guidewire in Khelix Sterilisation Tray A-T02 (or similar) **[image 8]**. Secure the handle/guidewire in the support brackets, ensuring the guidewire does not kink. Handle/guidewire assemblies are limited to one per tray.
2. Check that the reusable handle/guidewire is properly fitted into the support brackets and is not damaging the device.
3. Place tray lid on and lock the tray.
4. Wrap the tray in two (2) layers of 1-ply polypropylene wrap (e.g. Halyard Health H600 – 510(k) K082554 or similar) using sequential envelope folding techniques.
5. Placed the wrapped tray in the steriliser and sterilised using the validated washing parameters below:

Table 5: Validated Steam Sterilisation Parameters	
Steriliser Type	Prevacuum
Preconditioning Pulses	4
Temperature	132°C
Full Cycle Time	4 minutes
Dry Time	30 minutes
Device Configuration	Tray filled with one reusable handle and wrapped in two (2) layers of 1-ply polypropylene wrap (Halyard Health H600 –

	510(k) K082554) using sequential envelope folding techniques
Load Configuration	One 20-pound tray filled with miscellaneous stainless-steel dunnage placed above the tray containing the reusable handle in the steriliser.

6. Remove the sterilised wrapped tray from the steriliser and store according to storage instructions below.

The reprocessing instructions provided have been validated by CathRx Ltd. It remains the responsibility of the user to ensure reprocessing of the device is performed by trained personnel using the equipment and materials of the healthcare facility. Validation and routine monitoring of the process may be required to ensure desired results are achieved.

ASSEMBLY INSTRUCTIONS

The complete assembly procedure must be performed aseptically. No part of the sterile sheath or the sterile handle/guidewire assembly must contact any object outside the sterile field.

1. Match the handle/guidewire to the required catheter sheath (see **Table 2**). The deflection knob **[Image 9]** indicates the curve type **(I)** and the pin connector **(J)** of the handle/guidewire.
2. Before assembly, prepare an adequate sterile flat surface (e.g. sterile draped table) in the sterile field for the sheath and handle/guidewire.
3. Assemble immediately prior to use.
4. Remove the handle/guidewire from its sterile packaging or the sterilised wrapped tray using standard aseptic technique.
5. Check that the guidewire is fully un-deflected and straight, and ensure it is maintained in this position throughout the assembly process.

Always ensure the deflection knob is in the neutral position (sheath un-deflected) during assembly and disassembly to avoid damage to the guidewire.

6. Place the handle/guidewire on the flat sterile surface, with the guidewire coiled in one large loop **[Image 10]**.
7. Remove a new sheath from its sterile packaging and lay it straight on the flat sterile surface.
8. Hold handle and distal end of guidewire in one hand with the index finger and thumb approximately 2 cm from the guidewire tip **[Image 11]**.
9. Hold the sheath electrical connector with the other hand. Insert guidewire into centre hole of connector

and gently uncoil the guidewire to smoothly slide it all the way into the sheath [Image 12].

10. Insert the sheath electrical connector into the distal end of the handle and gently rotate until it is fully inserted and has clicked firmly into place [Image 13].
11. Check the catheter shaft can be smoothly deflected and straightened. Ensure that the catheter shaft is not coiled when this check is performed.

LIMITATIONS ON REPROCESSING

Reprocessing of Khelix Steerable Electrophysiology Catheter with Reuse Functionality has been validated for 10 cycles. One cycle includes cleaning and sterilisation processes.

STORAGE

Store in a cool, dry place (minimum temperature: -10°C [14°F]; maximum temperature: 55°C [131°F]; relative humidity: 85% ± 5%).

CONTAINMENT AND TRANSPORT

1. Remove the used catheter sheath and dispose. Refer to disassembly instructions provided above.
2. Wipe the used handle/guidewire with a lint-free cloth dampened with tap water to contain contaminants and prevent body fluids from drying.
3. Place and secure the used handle/guidewire inside Khelix Sterilisation Tray A-T02 (or similar). Place tray lid on and lock the tray.
4. Place the tray into an enclosed transport bin or enclosed case cart to prevent cross-contamination during transport.

DISPOSAL

Dispose device according to standard hospital procedure for biologically contaminated material.

WARNINGS AND PRECAUTIONS

- Do not attempt to operate the catheter prior to completely reading and understanding these Instructions for Use.
- Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement and placement should be done under fluoroscopic guidance.
- Do not use excessive force to advance or withdraw the catheter when resistance is encountered.
- The catheter should only be used by clinicians trained in cardiac electrophysiology procedures, in a fully equipped electrophysiology laboratory.
- The catheter sheath is for single use only. Do not reuse, reprocess or resterilise the sheath. Reuse, reprocessing or resterilisation may compromise the structural integrity of the device and/or lead to device failure. For example, detachment of small

components like tip, ring electrode dislocates and unable to sense signals in heart. Risk of cross contamination and poor catheter handling.

- Use of anti-coagulants should be considered for left-side procedures to reduce risk of thrombi formation and stroke.
- Cardiac catheterisation procedures present the potential for significant X-ray exposure. Exposure to X-rays can result in acute radiation injury as well as increased risk for somatic and genetic effects in both patients and electrophysiology laboratory staff, due to the X-ray beam intensity and duration of the fluoroscopic imaging. Steps should be taken to minimise X-ray exposure.
- In view of the potential for X-ray exposure and the associated risks, careful consideration must be given for the use of the device in pregnant women.
- Care should be taken when introducing, manipulating and withdrawing the catheter to avoid cardiac trauma and possible damage to the catheter which may result in detachment of catheter components.
- The catheter electrodes, the catheter and extension cable connectors, and the dispersive electrode, should not contact other conductive parts and earth.
- To avoid char formation on the catheter electrodes, do not apply RF energy when an ablation catheter is in contact with one or more of the Khelix Steerable Electrophysiology Catheter electrodes.
- Do not expose the catheter to organic solvents such as alcohol.
- When connected to power source, do not immerse the catheter handle or any cable connector in fluid as electrical performance could be affected.
- The catheter is MR unsafe. Do not use the catheter system near MRI equipment, as movement or heating of the catheter may occur, and the image or display may become distorted.
- The catheter is a type CF applied part. Any equipment used in connection with this device and its accessories must be defibrillation-proof type CF and meet IEC 60601-1 and IEC 60601-1-2.
- The catheter, together with its accessories, has been tested to allow connection with stimulation equipment with a rated current of 25mA, and a rated voltage of 27V.
- If the catheter becomes kinked, dispose of it. Replace with an undamaged device.
- Serious incidents related to the device should be reported to the manufacturer and the competent authority of your country.

WARRANTY AND DISCLAIMER

CathRx warrants that if it determines that a product was defective or faulty in materials and/or workmanship at

the time of purchase, CathRx will, at its option, provide a replacement CathRx product or refund the purchase price of the defective or faulty product.

This limited warranty only applies if the following conditions are met:

- The product was packaged and labelled by CathRx;
- The product is returned to CathRx for evaluation within 30 days of identification of the defect or fault by the original purchaser;
- The product has not been repaired, altered, modified, mishandled or reprocessed inappropriately;
- The product has been used, stored, cleaned, sterilised and reprocessed in accordance with the product labelling and these Instructions for Use; and
- The product is not used after the 'Use By' date marked on the packaging.

OTHER THAN AS EXPRESSLY SET OUT IN THIS LIMITED WARRANTY, TO THE MAXIMUM EXTENT PERMITTED BY LAW, CATHRX MAKES NO WARRANTY IN RELATION TO THE PRODUCT, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION ANY WARRANTY OF SAFETY, MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR THAT THE PRODUCT WILL BE DEFECT OR FAULT FREE. TO THE MAXIMUM EXTENT PERMITTED BY LAW, CATHRX EXCLUDES ALL SUCH WARRANTIES.

CathRx accepts no liability and the purchaser of the product assumes all liability, whether based on warranty, tort, contract, negligence, under statute or otherwise, for any kind of loss or damage (including special, incidental or consequential) or injury (including death) arising directly or indirectly as a result of or in connection with the products, including arising as a result of the handling, possession, use or misuse of the product. CathRx's entire liability, including for a breach of warranty or condition which cannot lawfully be modified or excluded, is limited, at the option of CathRx, to either the replacement of the CathRx product or the refund of the purchase price of the product. The remedies set forth in this limited warranty are the exclusive remedy available to any person. The use of the product is deemed to be acceptance of the terms and conditions of this limited warranty.

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