Diagnostic Extension Cable

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
Advena Ltd
Tower Business Centre
2nd Flr.
Tower Street
Swatar, BKR 4013
Malta

www.cathrx.com

CONTENTS

Images	2
English	3
Deutsch (German)	6
Français (French)	10
Italiano (Italian)	14
Nederlands (Dutch)	18
Español (Spanish)	22
Português (Portuguese)	26
Svenska (Swedish)	30
Čeština (Czech)	34
Dansk (Danish)	38
Norsk (Norwegian)	41
Magyar (Hungarian)	44
Symbol Definition	48

KHELIX DIAGNOSTIC EXTENSION CABLE

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- STERILE. Sterilized with Ethylene Oxide.
- Do not use if the package is open or damaged.
- Khelix Diagnostic Extension Cable is reusable per reprocessing instructions provided below.
- Use the device prior to the "Use By" date on the package label.

These instructions apply to the following products:

Product Name	Product Reference Number (REF)
Diagnostic Extension Cable	X-D10 X-D20

Sterilization: This product and its packaging have been sterilized with ethylene oxide gas (EO). Even though the product is processed in compliance with all applicable laws and regulations relating to EO exposure, Proposition 65, a State of California voter initiative, requires the following notice:

Warning: This product and its packaging have been sterilized with ethylene oxide. The packaging may expose you to ethylene oxide, a chemical known to the State of California to cause cancer or birth defects or other reproductive harm.

DEVICE DESCRIPTION

Khelix Diagnostic Extension Cable is an accessory for Khelix Diagnostic Electrophysiology Catheter.

The Khelix Diagnostic Extension Cable provides the interface between a Khelix Diagnostic Electrophysiology Catheter and standard electrophysiological recording equipment, such as ECG monitoring equipment, impedance based navigational equipment and cardiac stimulation equipment.

Khelix Diagnostic Extension Cable are available in both 10-pin for electrophysiology catheters with 4 to 10 electrodes and 20-pin for electrophysiology catheters with more than 10 electrodes. An example of 20-pin Khelix Diagnostic Extension Cable is shown in **[Image 1]**.

Khelix Diagnostic Extension Cable is intended to be used with the following devices:

Khelix Diagnostic Electrophysiology Catheters:
<ul style="list-style-type: none"> • Khelix Fixed Electrophysiology Catheter • Khelix Steerable Electrophysiology Catheter

Table 2: Devices intended to be used with Khelix Diagnostic Extension Cable

- Khelix Loop Electrophysiology Catheter
- Khelix Steerable Electrophysiology Catheter with Reuse Functionality

Khelix Sterilisation Tray (Ref: A-T01)

These devices are available and sold separately by CathRx Ltd.

INDICATIONS

Khelix Diagnostic Extension Cable provides interface between electrophysiology catheters and standard electrophysiological recording equipment, such as ECG monitoring equipment, impedance based navigational equipment and cardiac stimulation equipment.

INTENDED USER

Electrophysiologist, cardiac electrophysiologist, arrhythmia specialist or EP doctor with a specialization in abnormal heart rhythms. Electrophysiologists test the electrical activity of the heart to diagnose the source of arrhythmias (irregular heartbeats) to help determine a suitable treatment.

CONTRAINDICATIONS

There are no known contraindications for Diagnostic Extension Cable.

DIRECTIONS FOR USE

At point of use

1. Inspect the extension cable and packaging before opening. The contents of the package are sterile unless the package is opened or damaged.
2. If the package is opened or damaged, or if the package was opened and the extension cable unused, do not use the extension cable.

Start of Procedure

1. Remove new extension cable from its sterile packaging or the re-sterilised extension cable from the wrapped Khelix Sterilisation Tray A-T01 using aseptic technique. Place device on a sterile working area.
2. Inspect the extension cable for any damages including damage to the insulation or connectors. The connector contacts must be dry.
3. Connect the extension cable to the catheter by inserting the distal end of the cable into the socket located at the proximal end of the catheter handle. Prior to connection, ensure the key (raised ridges) on the cable connector (**A**) are aligned with the key slots inside the connector socket of the catheter handle as shown in **[Image 2]**. To prevent catheter shaft deflection during connection, grip the catheter

deflection knob **(B)** during cable insertion as shown in **[Image 3]**.

4. Connect the extension cable to the system.

End of Procedure

1. Disconnect extension cable from catheter by gripping the extension cable pin head **(C)** where the arrows are printed on the pin head and pull outwards as shown in **[Image 2]**. For deflectable catheters, to prevent catheter shaft deflection during disconnection, grip the catheter deflection knob **(B)** during disconnection as shown in **[Image 3]**.
2. Disconnect extension cable from system, immediately wipe down the cable with a sponge or towel moistened with water.
3. If not reusing, dispose the extension cable according to standard hospital procedure for biologically contaminated material.

REPROCESSING INSTRUCTIONS

-  The extension cable can be reprocessed following the instructions below, provided that the limitations for reprocessing indicated below is not exceeded.

Manual Cleaning Method

1. Prepare ENZOL® Enzymatic Detergent solution (or equivalent) in a container by diluting 30 ml of detergent in 3785 ml of tap water (20–25°C).

 Wear gloves during solution preparation and cleaning procedures.

2. Rinse the extension cable under running tap water (20–25°C) for 30 seconds.
3. Soak the extension cable in the detergent solution for 5 minutes.
4. Remove the extension cable from the detergent solution and place on a clean CSR (sterilisation) drape.
5. Using a non-abrasive nylon bristled brush saturated with the detergent solution, gently brush all surfaces of the device for 2 minutes.
6. Rinse device in running tap water (20–25°C) for 30 seconds.
7. Immerse device in de-ionized (DI) water 3 separate times for 10 seconds each.
8. Wipe the device using a non-linting wipe. Place on a clean CSR (sterilisation) drape and allow to dry until visibly free of moisture.
9. After drying, inspect visually for evidence of soil remaining on device. If visible soil is present, remove soil by manual cleaning and brushing using soft sponge/cloth or soft-bristled brush, and/or repeat of manual cleaning prior to sending to sterilisation.

Automatic Cleaning Method

1. Place the extension cable in the washer equipment manufacturer's tray or a suitable CSD washing tray.
2. Clean the extension cable in a Bellimed Washer Model WD250 using Ecolab® Neutral Enzymatic Detergent (or an equivalent automated washer/enzymatic detergent) with the following validated washing parameters:
 - Wash at 45°C for 4 minutes. Dose pump 4 (detergent) 5 ml
 - Wash at 60°C for 3 minutes
 - Rinse with unheated water for 1 minute
 - Rinse at 60°C for 1 minute
 - Thermal Disinfection A₀ 93°C 2.5 minutes
 - A₀ value: A₀3000
 - Dry at 123°C air for 14 minutes.
3. Upon completion of wash cycle, use gloves to remove the device. Wipe using a non-linting wipe. Place on a clean CSR (sterilisation) drape. Dry until visibly free of moisture.
4. After drying, inspect visually for evidence of soil remaining on each device. If visible soil is present, remove soil by manual cleaning and brushing using soft sponge/cloth or soft-bristled brush, and/or repeat of automatic cycles prior to sending to sterilisation.

Failure to properly clean the device could lead to inadequate sterilisation.

Inspection and Maintenance

1. Visually inspect the extension cable for any damages.
2. If extension cable is not damaged or kinked, continue the use and proceed to sterilisation.
3. Dispose any damaged extension cable according to disposal instructions provided below.

Steam Sterilisation

1. Place extension cable in Khelix Sterilisation Tray A-T01 **[Image 8]**. Secure extension cable in holders as indicated by the markings within the tray. Extension cable is limited to one per tray.
2. Wrap the trays in a sterilisation wrap (eg Kimberly-Clark Kinguard One-Step Sterilization Wraps) using a simultaneous double-wrapping envelope fold technique.
3. The following steam sterilisation parameters were validated in an AMSCO Vacumatic "A" Autoclave:
 - Cycle Type: Dynamic Air Removal
 - Pre-vacuum pressure: 27" Hg (0.914 bar)
 - Pre-vacuum Pulses: 3
 - Dwell Temperature: 132°C
 - Dwell Pressure: 26.8 psig (1.85 bar)
 - Dwell Time: 4 minutes
 - Drying Time: 20 minutes.

The reprocessing instructions provided have been validated by CathRx Ltd. It remains the responsibility of the user to ensure reprocessing of the device is performed by trained personnel using the equipment and materials of the healthcare facility. Validation and routine monitoring of the process may be required to ensure desired results are achieved.

LIMITATIONS ON REPROCESSING

Reprocessing Khelix Diagnostic Extension Cable has been validated for 20 cycles. One cycle consists of cleaning and sterilisation processes.

STORAGE

Store in a cool, dry place (minimum temperature: -10°C [14°F]; maximum temperature: 55°C [131°F]; relative humidity: 85% ± 5%).

CONTAINMENT AND TRANSPORT

1. Disconnect the extension cable from the system, wipe the extension cable with a lint-free cloth dampened with tap water to contain contaminants and prevent body fluids from drying.
2. Place and secure the used extension cable inside Khelix Sterilisation Tray A-T01. Place tray lid on and lock the tray.
3. Place the tray into an enclosed transport bin or enclosed case cart to prevent cross-contamination during transport.

DISPOSAL

Dispose device according to standard hospital procedure for biologically contaminated material.

WARNINGS AND PRECAUTIONS

- The extension cable connectors should not contact other conductive parts and earth.
- Do not expose the extension cable to organic solvents such as alcohol.
- When connected to power source, do not immerse any cable connector in fluid as electrical performance could be affected.
- The extension cable is a type CF applied part. Any equipment used in connection with this device and its accessories must be defibrillation-proof type CF and meet IEC 60601-1 and IEC 60601-1-2.
- The extension cable, together with the catheters, has been tested to allow connection with stimulation equipment with a rated current of 25mA, and a rated voltage of 27V.
- If the extension cable is damaged, dispose of it. Replace with an undamaged device.

- Serious incidents related to the device should be reported to the manufacturer and the competent authority of your country.

WARRANTY AND DISCLAIMER

CathRx warrants that if it determines that a product was defective or faulty in materials and/or workmanship at the time of purchase, CathRx will, at its option, provide a replacement CathRx product or refund the purchase price of the defective or faulty product.

This limited warranty only applies if the following conditions are met:

- The product was packaged and labelled by CathRx;
- The product is returned to CathRx for evaluation within 30 days of identification of the defect or fault by the original purchaser;
- The product has not been repaired, altered, modified, mishandled or reprocessed inappropriately;
- The product has been used, stored, cleaned, sterilised and reprocessed in accordance with the product labelling and these Instructions for Use; and
- The product is not used after the 'Use By' date marked on the packaging.

OTHER THAN AS EXPRESSLY SET OUT IN THIS LIMITED WARRANTY, TO THE MAXIMUM EXTENT PERMITTED BY LAW, CATHRX MAKES NO WARRANTY IN RELATION TO THE PRODUCT, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION ANY WARRANTY OF SAFETY, MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR THAT THE PRODUCT WILL BE DEFECT OR FAULT FREE. TO THE MAXIMUM EXTENT PERMITTED BY LAW, CATHRX EXCLUDES ALL SUCH WARRANTIES.

CathRx accepts no liability and the purchaser of the product assumes all liability, whether based on warranty, tort, contract, negligence, under statute or otherwise, for any kind of loss or damage (including special, incidental or consequential) or injury (including death) arising directly or indirectly as a result of or in connection with the products, including arising as a result of the handling, possession, use or misuse of the product. CathRx's entire liability, including for a breach of warranty or condition which cannot lawfully be modified or excluded, is limited, at the option of CathRx, to either the replacement of the CathRx product or the refund of the purchase price of the product. The remedies set forth in this limited warranty are the exclusive remedy available to any person. The use of the product is deemed to be acceptance of the terms and conditions of this limited warranty.

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