# IFU - Khelix Fixed Electrophysiology Catheter

Doc # IFU0008  
Rev: V  
Status: Live

Revision	Originator	ECR #	Release date	Summary of Change	Reason for Change
Q	L.Huang	ECR00589	11 Jul 2018	Update IFU to add extension cable connection and disconnection instruction in English and German for clarification at beginning and end of the procedure.	Emphasize cable/catheter pins to be aligned during extension cable connection. Cable pin head to be gripped and pulled away for easy cable disconnection after procedure.
R	L.Huang	ECR00641	07 Jan 2019	Reformat IFU0008 to add 9 more additional languages. French, Italian, Dutch, Spanish, Portuguese, Swedish, Czech, Danish, Norwegian. Remove non-mandatory symbol definition. Add latex free symbol.	Reformat IFU0008 to include more languages. To define only mandatory symbols required to be explained in information supplied by manufacturer. Marketing purpose to be in line with competitor's product.
S	L.Huang	ECR00694	1 Apr 2019	Change EU authorised representative from Priority Analyst Ltd to Advena Ltd whom is based in Malta. Update NB0086 to NB2797.	UK leaving EU, Brexit.
T	C. Cuizon	ECR00800	27 Feb 2020	1. Reformatted. 2. Updated information supplied per applicable standards and regulatory requirements.	1. Clearer flow of information. Good documentation practice. 2. Comply to MDR and FDA labelling requirements and other applicable standards (ISO 17664:2017, etc.)
U	C. Cuizon	ECR00842	3 Aug 2020	Updated in Table 1 to include prefix of applicable devices.	Updated Product reference column (Table 1) to reflect the general prefix of all applicable devices. This IFU applies to Fixed EP Catheter range with product reference prefix [5SSF-*]. Note: No change to the contents e.g. indication, contraindications, warnings/precautions and directions for use, etc.
V	I. Alkhalaeileh	ECR00920	20 Jan 2022	1. Modified storage conditions to specify only "Keep dry. Keep away from sunlight." for all languages. 2. Minor formatting changes. 3. Changed date in footer to an automatic field.	1. Changes were made to clarify storage conditions for device. In addition, CathRx devices have been tested to this specification, however it is not necessary for the storage conditions. 2. Changes were made due to the minor size changes the modification caused. 3. This change was made to ensure that the release date was correctly maintained.

# Fixed Electrophysiology Catheter

Fixierte Elektrophysiologie-Katheter | Cathéter Électrophysiologique Fixe | Catetere Fisso Per Elettrofisiologia | Gefixeerde Elektrofysiologiekatheter | Catéter de Electrofisiologia Fijo | Cateter de Eletrofisiologia Fijo | Fasta Electrofysiologikateter | Fixni Elektrofiziologický Katetr | Fast Elektrofysiologisk Kateter | Fiksert Elektrofysiologikateter

## Instructions for use

Gebruiksaanweisung | Conseils d'utilisation | Istruzioni d'uso | Gebruiksaanwijzingen | Instrucciones de uso | Instruções de uso | Bruksanvisning | Návod k použití | Brugsanvisning | Bruksanvisning

CathRx Ltd  
Unit 8/2-8 South Street  
Rydalmere NSW 2116  
Australia  
Tel: (+61) 2 7804 1846  
Email: [info@cathrx.com](mailto:info@cathrx.com)  
[www.cathrx.com](http://www.cathrx.com)

Advena Ltd  
Tower Business Centre  
2<sup>nd</sup> Flr, Tower Street  
Swatar, BKR 4013  
Malta

[www.cathrx.com](http://www.cathrx.com)

## CONTENTS

English	2
Deutsch (German)	5
Français (French)	8
Italiano (Italian)	11
Nederlands (Dutch)	14
Español (Spanish)	17
Português (Portuguese)	20
Svenska (Swedish)	23
Čeština (Czech)	26
Dansk (Danish)	29
Norsk (Norwegian)	32
Symbol Definition	35

## KHELIX FIXED ELECTROPHYSIOLOGY CATHETER

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- STERILE. Sterilized with Ethylene Oxide.
- Do not use if the package is open or damaged.
- The catheter is for single use only. Do not reuse, reprocess or re-sterilise the devices.
- Use the device prior to the "Use By" date on the package label.

These instructions apply to the following products:

Table 1: Product Names

Product Name	Product Reference Number (REF)
Fixed Electrophysiology Catheter	Ref: [5SSF-*.*]

**Sterilization:** This product and its packaging have been sterilized with ethylene oxide gas (EO). Even though the product is processed in compliance with all applicable laws and regulations relating to EO exposure, Proposition 65, a State of California voter initiative, requires the following notice:

Warning: This product and its packaging have been sterilized with ethylene oxide. The packaging may expose you to ethylene oxide, a chemical known to the State of California to cause cancer or birth defects or other reproductive harm.

### DEVICE DESCRIPTION

The Khelix Fixed Electrophysiology Catheter is designed to facilitate electrophysiological mapping of the heart. It consists of a catheter shaft containing an array of platinum-iridium electrodes for stimulation and recording. The device is available in a number of different electrode configurations and curve types.

The Khelix Fixed Electrophysiology Catheter interfaces with standard electrophysiological recording equipment, such as ECG monitoring equipment, impedance based navigational equipment and cardiac stimulation equipment, via the Khelix Diagnostic Extension Cable (supplied separately).

Khelix Fixed Electrophysiology Catheter is intended to be used with the following accessory:

Table 2: Accessory for Khelix Fixed Electrophysiology Catheter
Khelix Diagnostic Extension Cable

This accessory is available and sold separately by CathRx Ltd.

### INDICATIONS

The Khelix Fixed Electrophysiology Catheter is indicated for the electrical recording or stimulation of endocardial structures.

### INTENDED USER

Khelix Fixed Electrophysiology Catheter is used by interventional cardiologists in electrophysiology (EP) studies to determine the cause of an abnormal heart rhythm, to locate the site of origin of an abnormal heart rhythm, to decide the best treatment for an abnormal heart rhythm, and/or to check the effectiveness of ablation therapy.

### CONTRAINDICATIONS

- The catheter has not been shown to be safe and effective for radio frequency (RF) ablation.
- The transeptal approach is contraindicated in patients with left atrial thrombus, myxoma, or an inter-atrial patch or baffle.
- The use of this catheter may not be appropriate for use in patients with prosthetic valves.
- Electrophysiology procedures are contraindicated for patients with unstable cardiac conditions, e.g. acute myocardial infarction, unstable angina, hemodynamic instability.
- There is a relative contraindication for cardiac catheterisation procedures in patients with active systemic infection.

### ADVERSE REACTIONS

A number of serious adverse reactions have been documented for cardiac catheterisation procedures, including cardiac tamponade, pulmonary embolism, myocardial infarction, stroke and death.

The following additional complications associated with cardiac catheterisation have also been reported in the literature: vascular bleeding, local haematoma, thrombosis, AV fistula, pseudoaneurysm, thromboembolism, vasovagal reaction, cardiac perforation, air embolism, arrhythmia, valvular damage, pneumothorax and haemothorax.

### DIRECTIONS FOR USE

#### At point of use

1. Inspect the catheter and packaging before opening. The contents of the package are sterile unless the package is opened or damaged.
2. If the package is opened or damaged, or if the package was opened and the catheter unused, do not use the catheter. Do not attempt to resterilise.
3. A detachable package label has been provided and may be affixed to the medical record of the patient.

**Start of Procedure**

1. Remove the catheter from its sterile packaging using aseptic technique. Place device on a sterile working area.
2. Inspect the catheter for any damages. Do not use the catheter if kinked.
3. Inspect the catheter for damage to the insulation or connectors. The connector contacts must be dry.
4. Test all functions of the catheter.
5. Connect the extension cable by inserting the distal end of the extension cable into the socket located at the proximal end of the catheter. Prior to connection, ensure the key (raised ridges) on the cable connector are aligned with the key slots inside the connector socket of the catheter.

The catheter must only be used with Khelix Diagnostic Extension Cable.

**During Procedure**

1. Vascular access is created with a minimum 6 French haemostasis introducer sheath to accommodate the catheter.
2. Advance the catheter through the vascular access site into the desired intracardiac position using fluoroscopic and ECG guidance.

**End of Procedure**

1. Remove the catheter from the patient.
2. Disconnect extension cable from catheter by gripping the extension cable pin head where the arrows are printed on the pin head and pull outwards.
3. Dispose the catheter according to standard hospital procedure for biologically contaminated material.

**STORAGE**

Keep dry.

Keep away from sunlight.

**DISPOSAL**

Dispose device according to standard hospital procedure for biologically contaminated material.

**WARNINGS AND PRECAUTIONS**

- Do not attempt to operate the catheter prior to completely reading and understanding these Instructions for Use.
- Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement and placement should be done under fluoroscopic guidance.
- Do not use excessive force to advance or withdraw the catheter when resistance is encountered.
- The catheter should only be used by clinicians trained in cardiac electrophysiology procedures, in a fully equipped electrophysiology laboratory.
- The catheter is for single use only. Do not reuse, reprocess or resterilise the device. Reuse, reprocessing or resterilisation may compromise the structural integrity of the devices and/or lead to device failure. For example, detachment of small components like tip, ring electrode dislocates, unable to sense signals in heart. Risk of cross contamination and poor catheter handling.
- Use of anti-coagulants should be considered for left-side procedures to reduce risk of thrombi formation and stroke.
- Cardiac catheterisation procedures present the potential for significant X-ray exposure. Exposure to X-rays can result in acute radiation injury as well as increased risk for somatic and genetic effects in both patients and electrophysiology laboratory staff, due to the X-ray beam intensity and duration of the fluoroscopic imaging. Steps should be taken to minimise X-ray exposure.
- In view of the potential for X-ray exposure and the associated risks, careful consideration must be given for the use of the device in pregnant women.
- Care should be taken when introducing, manipulating and withdrawing the catheter to avoid cardiac trauma and possible damage to the catheter which may result in detachment of catheter components.
- The catheter electrodes, the catheter and extension cable connectors, and the dispersive electrode, should not contact other conductive parts and earth.
- To avoid char formation on the Khelix Fixed Electrophysiology Catheter electrodes, do not apply RF energy when an ablation catheter is in contact with one or more of the catheter electrodes.
- Do not expose the catheter to organic solvents such as alcohol.
- Do not immerse the catheter handle or any cable connector in fluid as electrical performance could be affected.
- The catheter is MR unsafe. Do not use the catheter system near MRI equipment, as movement or heating of the catheter may occur, and the image or display may become distorted.
- The catheter is a type CF applied part. Any equipment used in connection with this device and its accessories must be defibrillation-proof type CF and meet IEC 60601-1 and IEC 60601-1-2.
- The catheter, together with its accessories, has been tested to allow connection with stimulation equipment with a rated current of 25mA, and a rated voltage of 27V.
- If the catheter becomes kinked, dispose of it. Replace with an undamaged device.
- Serious incidents related to the device should be reported to the manufacturer and the competent authority of your country.

## **WARRANTY AND DISCLAIMER**

CathRx warrants that if it determines that a product was defective or faulty in materials and/or workmanship at the time of purchase, CathRx will, at its option, provide a replacement CathRx product or refund the purchase price of the defective or faulty product.

This limited warranty only applies if the following conditions are met:

- The product was packaged and labelled by CathRx;
- The product is returned to CathRx for evaluation within 30 days of identification of the defect or fault by the original purchaser;
- The product has not been repaired, altered, modified, mishandled or reprocessed inappropriately;
- The product has been used, stored, cleaned, sterilised and reprocessed in accordance with the product labelling and these Instructions for Use; and
- The product is not used after the 'Use By' date marked on the packaging.

OTHER THAN AS EXPRESSLY SET OUT IN THIS LIMITED WARRANTY, TO THE MAXIMUM EXTENT PERMITTED BY LAW, CATHRX MAKES NO WARRANTY IN RELATION TO THE PRODUCT, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION ANY WARRANTY OF SAFETY, MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR THAT THE PRODUCT WILL BE DEFECT OR FAULT FREE. TO THE MAXIMUM EXTENT PERMITTED BY LAW, CATHRX EXCLUDES ALL SUCH WARRANTIES.

CathRx accepts no liability and the purchaser of the product assumes all liability, whether based on warranty, tort, contract, negligence, under statute or otherwise, for any kind of loss or damage (including special, incidental or consequential) or injury (including death) arising directly or indirectly as a result of or in connection with the products, including arising as a result of the handling, possession, use or misuse of the product. CathRx's entire liability, including for a breach of warranty or condition which cannot lawfully be modified or excluded, is limited, at the option of CathRx, to either the replacement of the CathRx product or the refund of the purchase price of the product. The remedies set forth in this limited warranty are the exclusive remedy available to any person. The use of the product is deemed to be acceptance of the terms and conditions of this limited warranty.

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