Sterilisation Tray for Reusable Handle

Sterilisationstablett für wiederverwendbaren Handgriff | Bac de stérilisation pour poignée réutilisable | Vassoio di sterilizzazione per l'impugnatura riutilizzabile | Sterilisatiebak voor herbruikbare handgreep | Bandeja de esterilización para mango reutilizable | Bandeja de esterilização para manopla reutilizável | Steriliseringskorg för återanvändbart handtag | Sterilizační sito na rukojeť k opakovanému použití | Steriliseringsbakke til genanvendelige håndtag | Steriliseringsbrett for gjenbrukbart håndtak

Instructions for use

Gebrauchsanweisung | Conseils d'utilisation | Istruzioni d'uso | Gebruiksaanwijzingen | Instrucciones de uso | Instruções de uso | Bruksanvisning | Návod k použití | Brugsanvisning | Bruksanvisning

CathRx Ltd
Unit 8/2-8 South Street
Rydalmere NSW 2116
Australia
Tel: (+61) 2 7804 1846
Email: info@cathrx.com
www.cathrx.com

Advena Ltd
Tower Business Centre
2nd Flr, Tower Street
Swatar, BKR 4013
Malta

www.cathrx.com

CONTENTS

| | |
|------------------------|----|
| Images | 2 |
| English | 3 |
| Deutsch (German) | 6 |
| Français (French) | 10 |
| Italiano (Italian) | 14 |
| Nederlands (Dutch) | 18 |
| Español (Spanish) | 22 |
| Português (Portuguese) | 26 |
| Svenska (Swedish) | 30 |
| Čeština (Czech) | 34 |
| Dansk (Danish) | 37 |
| Norsk (Norwegian) | 41 |
| Symbol Definition | 44 |

Image 1: Sterilisation Tray A-T02¹ (DIN1/1 Type Tray)



Image 2: Top View securing Khelix Steerable Electrophysiology Catheter with Reuse Functionality (Reusable Handle/Guidewire)

| | |
|--|-------------------------|
| A. Tray Handles (2) | B. Support brackets (7) |
| C. Khelix Steerable Electrophysiology Catheter with Reuse Functionality ² (handle/guidewire assembly) | |

¹ A-T02 Sterilisation Tray can also be ordered directly from supplier Van Straten Medical under order code: 262601639. Website: www.vanstratenmedical.com

² Khelix Steerable Electrophysiology Catheter with Reuse Functionality is supplied separately and available for purchase at CathRx Ltd.

KHELIX STERILISATION TRAY FOR REUSABLE HANDLE

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- Khelix Sterilisation Tray for Reusable Handle is supplied non-sterile and is intended for re-use per reprocessing instructions provided below.

These instructions apply to the following products:

| Table 1: Product Name | |
|--|--------------------------|
| Product Name | Product Reference Number |
| Sterilisation Tray for Reusable Handle | A-T02 |

DEVICE DESCRIPTION

Khelix Sterilisation Tray for Reusable Handle is a stainless-steel tray composed of a body and lid with perforated holes to facilitate access of steam during autoclave. The tray is a DIN 1/1 type tray with dimensions of 480x250x60mm (LxWxH) **[Image 1]**. Brackets allow securement of Khelix Steerable Electrophysiology Catheter with Reuse Functionality within the tray body **[Image 2]**. The support brackets are made up of silicone and the tray handle is made up of stainless-steel tray.

The tray is used in conjunction with a sterilisation wrap as a sterile barrier during and after steam sterilisation. Refer to sterilisation wrap manufacturer's instructions for use for shelf-life information.

Khelix Sterilisation Tray for Reusable Handle is intended to be used with the following devices:

| Table 2. Devices to be used with Khelix Sterilisation Tray for Reusable Handle |
|--|
| Khelix Steerable Electrophysiology Catheter with Reuse Functionality* |

*Also referred as reusable handle/guidewire in this IFU.

INDICATIONS

Khelix Sterilisation Tray for Reusable Handle is intended to secure and protect the reusable handle/guidewire assembly of Khelix Steerable Electrophysiology Catheter with Reuse Functionality (after removal of the catheter sheath. Refer to Khelix Steerable Electrophysiology Catheter with Reuse Functionality instructions for use for the disassembly procedure) during automated cleaning and steam sterilisation conducted by the hospital's central sterilisation department (CSD), and before and after procedure during device transportation between CSD and electrophysiology laboratory.

INTENDED USER

Khelix Sterilisation Tray for Reusable Handle is intended to be used by the hospital's central sterilisation department (CSD) during automated cleaning and steam sterilisation and transportation of the reusable handle/guidewire.

CONTRAINDICATIONS

There are no known contraindications for Khelix Sterilisation Tray for Reusable Handle.

DIRECTIONS FOR USE

At point of use

1. Check that tray is not damaged. Inspect tray components for cleanliness and confirm no visual contamination prior to use.
2. Refer to Khelix Steerable Electrophysiology Catheter with Reuse Functionality instructions for use for more information on handling the reusable handle/guidewire during reprocessing.
3. Place the reusable handle/guidewire inside the tray and secure the handle/guidewire in the support bracket as indicated by **[Image 2]**.
4. Clean, sterilise and transport sterilisation tray according to instructions provided below.

REPROCESSING INSTRUCTIONS

Refer to Khelix Steerable Electrophysiology Catheter with Reuse Functionality instructions for use for more information on re-processing of reusable handle/guidewire.

Automatic Cleaning Method

Note: Automated cleaning can be done for tray alone or with reusable handle/guidewire within the tray. If tray is to be cleaned alone, proceed to step 3.

1. Place the reusable handle/guidewire in the Khelix Sterilisation Tray (or similar) **[Image 2]**. Secure the reusable handle/guidewire in the support brackets, ensuring the guidewire does not kink. Handle/guidewire assemblies are limited to one per tray.
2. Check that the reusable handle/guidewire is properly fitted into the support brackets and is not damaging the device.
3. Place tray lid on and lock the tray.
4. Place the tray in a medical grade washer (e.g. Steris AMSCO models or Belimed WD models).
5. Start the automated cleaning using the validated washing parameters below:

| Table 3. Validated Automated Cleaning Parameters | | |
|--|----------------|-----------------------------|
| Wash 1 | Detergent Type | Neutral enzymatic detergent |
| | Water Quality | Tap Water |

Table 3. Validated Automated Cleaning Parameters

| | | |
|----------|--------------------|-----------------------------|
| | Volume (dose pump) | ¼ oz per gallon (1.87 g/l) |
| | Temperature | 45°C tap water (set point) |
| | Minimum Time | 4 minutes |
| Wash 2 | Detergent Type | Neutral enzymatic detergent |
| | Water Quality | Tap Water |
| | Volume (dose pump) | ¼ oz per gallon (1.87 g/l) |
| | Temperature | 60°C tap water (set point) |
| | Minimum Time | 3 minutes |
| Rinse | Water Quality | Tap Water |
| | Temperature | 60°C tap water (set point) |
| | Minimum Time | 1 minute |
| Dry Time | Temperature | 115°C (set point) |
| | Minimum Time | 14 minutes |

6. After completion of automated cleaning, visually inspect the tray for any remaining visible soil and contaminants. If visible soil is present, remove soil by manual cleaning and brushing using soft sponge/cloth or soft-bristled brush, and/or repeat of automatic cycles prior to sending to sterilisation.

Failure to properly clean the device could lead to inadequate sterilisation.

Visual Inspection

1. Visually inspect the tray for any damage. Ensure the lid, bottom and corners of the tray are not bent or dented. Ensure the locking mechanisms are in working order. Ensure the handles are not broken. Ensure the brackets are not broken or bent, ripped or displaced.
2. If no damage is observed, continue the use of tray and proceed to sterilisation.
3. Dispose any damaged tray per disposal instructions provided below and replace with new tray and repeat washing cycle for soiled catheter.

Steam Sterilisation

Note: skip step 1 to 3 below if the reusable handle/guidewire is already placed into the tray during automated washing cycle.

1. Place the cleaned reusable handle/guidewire in Khelix Sterilisation Tray (or similar) **[Image 2]**. Secure the reusable handle/guidewire in the support brackets, ensuring the guidewire does not kink. Handle/guidewire assemblies are limited to one per tray.
2. Check that the reusable handle/guidewire is properly fitted into the support brackets and is not damaging the device.

3. Place tray lid on and lock the tray.
4. Wrap the tray in two (2) layers of 1-ply polypropylene wrap (Halyard Health H600 – 510(k) K082554 or similar) using sequential envelope folding technique.
5. Placed the wrapped tray in the steriliser and sterilise according to parameters in Table 4.

Table 4. Validated Steam Sterilisation Parameters

| | |
|------------------------|--|
| Steriliser Type | Prevacuum |
| Preconditioning Pulses | 4 |
| Temperature | 132°C |
| Full Cycle Time | 4 minutes |
| Dry Time | 30 minutes |
| Device Configuration | Tray loaded with one reusable handle and wrapped in two (2) layers of 1-ply polypropylene wrap (Halyard Health H600 – 510(k) K082554) using sequential envelope folding techniques |
| Load Configuration | One 20-pound tray filled with miscellaneous stainless-steel dunnage placed above the tray containing the reusable handle in the steriliser. |

6. Remove the sterilised wrapped tray from the steriliser and store according to storage instructions below.

The reprocessing instructions provided have been validated by CathRx Ltd. It remains the responsibility of the user to ensure reprocessing of the device is performed by trained personnel using the equipment and materials of the healthcare facility. Validation and routine monitoring of the process may be required to ensure desired results are achieved.

Post-process handling

1. Prior to unwrapping the sterilised wrapped tray for clinical use, inspect the outside of the wrap to ensure the wrap is intact and has no tears or holes.

Unwrapping of sterilised tray should be done in a sterile field. Refer to sterilisation wrap manufacturer's instructions for use.
2. Check the tray for any damages.

LIMITATIONS ON REPROCESSING

The tray can be reused until unacceptable deterioration such as corrosion, cracking, rust, peeling, flaking, or mechanical failure occurs (see visual inspection criteria). Signs of mechanical failure include:

- Broken or cracked corners
- Broken handles
- Broken or ripped support brackets

STORAGE

- Keep dry.
- Keep away from sunlight.
- Refer to sterilisation wrap manufacturer's instructions for use for shelf life information.

CONTAINMENT AND TRANSPORT

Place and secure tray with or without Khelix Steerable Electrophysiology Catheter with Reuse Functionality into an enclosed transport bin or enclosed case cart to prevent cross-contamination during transport.

DISPOSAL

Dispose device according to standard hospital procedure for potential biologically contaminated material.

WARNINGS AND PRECAUTIONS

- The reprocessing instructions provided are validated specifically for Khelix Steerable Electrophysiology Catheter with Reuse Functionality and Khelix Sterilisation Tray for Reusable Handle.
- Do not overload tray. One tray is limited to one reusable handle/guidewire only.
- Do not use tray if it is damaged or broken. Dispose damaged or broken tray per disposal instructions below.
- Do not use abrasive cleaners or pads, metal or wire brushes or solvents when cleaning the tray.
- Khelix Sterilisation Tray for Reusable Handle is not designed to maintain sterility. They are designed to facilitate sterilisation process when used in conjunction with sterilisation wrap as a sterile barrier.
- Serious incidents related to the device should be reported to the manufacturer and the competent authority of your country.

WARRANTY AND DISCLAIMER

CathRx warrants that if it determines that a product was defective or faulty in materials and/or workmanship at the time of purchase, CathRx will, at its option, provide a replacement CathRx product or refund the purchase price of the defective or faulty product.

This limited warranty only applies if the following conditions are met:

- The product was packaged and labelled by CathRx;
- The product is returned to CathRx for evaluation within 30 days of identification of the defect or fault by the original purchaser;
- The product has not been repaired, altered, modified, mishandled or reprocessed inappropriately;

- The product has been used, stored, cleaned, sterilised and reprocessed in accordance with the product labelling and these Instructions for Use; and
- The product is not used after the 'Use By' date marked on the packaging.

OTHER THAN AS EXPRESSLY SET OUT IN THIS LIMITED WARRANTY, TO THE MAXIMUM EXTENT PERMITTED BY LAW, CATHRX MAKES NO WARRANTY IN RELATION TO THE PRODUCT, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION ANY WARRANTY OF SAFETY, MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR THAT THE PRODUCT WILL BE DEFECT OR FAULT FREE. TO THE MAXIMUM EXTENT PERMITTED BY LAW, CATHRX EXCLUDES ALL SUCH WARRANTIES.

CathRx accepts no liability and the purchaser of the product assumes all liability, whether based on warranty, tort, contract, negligence, under statute or otherwise, for any kind of loss or damage (including special, incidental or consequential) or injury (including death) arising directly or indirectly as a result of or in connection with the products, including arising as a result of the handling, possession, use or misuse of the product. CathRx's entire liability, including for a breach of warranty or condition which cannot lawfully be modified or excluded, is limited, at the option of CathRx, to either the replacement of the CathRx product or the refund of the purchase price of the product. The remedies set forth in this limited warranty are the exclusive remedy available to any person. The use of the product is deemed to be acceptance of the terms and conditions of this limited warranty.

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