

Khelix



European Catalogue

Table of Contents

Table of Contents	2
Product Types	3
Curve Types	4
Fixed Electrophysiology Catheter Curve Types	4
Steerable Electrophysiology Catheter Curve Types	9
Steerable Loop Curve Electrophysiology Catheter Types	11
Catheter Configurations	12
Khelix Fixed Curve Electrophysiology Catheters	12
Steerable Electrophysiology Catheters	18
Steerable Electrophysiology Catheter with Reuse Functionality	39
Disposable Sheaths for Catheters with Reuse Functionality	49
Loop Steerable Electrophysiology Catheter	51
Accessories	54
Sterilisation Trays	54
Extension Cables	54

Ordering Information

CathRx Ltd
Unit 8/2-8 South Street
Rydalmere NSW 2116
Australia
Tel: (+61) 2 7804 1846
Email: info@cathrx.com
www.cathrx.com

Product Types

FIXED CURVE
ELECTROPHYSIOLOGY
CATHETER



STEERABLE
ELECTROPHYSIOLOGY
CATHETER




STEERABLE
ELECTROPHYSIOLOGY
CATHETER




WITH REUSE
FUNCTIONALITY

LOOP STEERABLE
ELECTROPHYSIOLOGY
CATHETER




ACCESSORIES



Fixed Electrophysiology Catheter Types

Typical Fixed Curves - a comparison


The fixed curve catheter family is comprised of four curve types depicted below (typ), the following pages show each available curve variation.



Fixed Electrophysiology Catheter Types

A Curve Variations


CathRx Shape Name	Curve Size	Curve Radius R1	Length L1 (reach)
A Nominal reach (nom)	Nominal	32.0mm	0mm
A Extended reach	Nominal	32.0mm	10mm
A-Large Nominal reach	Large	38.0mm	0mm
A-Large Extended reach	Large	38.0mm	10mm



Fixed Electrophysiology Catheter Types

D Curve Variations


CathRx Shape Name	Curve Size	Curve Radius R1	Length L1 (reach)
D Nominal reach (nom)	Nominal	21.5mm	21mm
D Extended reach	Nominal	21.5mm	31mm
D-Large Nominal reach	Large	27.5mm	21mm
D-Large Extended reach	Large	27.5mm	31mm



Fixed Electrophysiology Catheter Types

F Curve Variations


CathRx Shape Name	Curve Size	Curve Radius R1	Length L1 (reach)
F Nominal reach (nom)	Nominal	29.5mm	15mm
F Extended reach	Nominal	29.5mm	25mm
F-Large Nominal reach	Large	35.5mm	15mm
F-Large Extended reach	Large	35.5mm	25mm



Fixed Electrophysiology Catheter Types

Y Curve Variations


CathRx Shape Name	Curve Size	Curve Radius R1	Length L1	Length L2
Y Nominal reach (nom)	Nominal	28.0mm	16mm	0mm
Y Extended reach	Nominal	28.0mm	26mm	6mm
Y X-Extended reach	Nominal	28.0mm	36mm	12mm
Y-Large Nominal reach	Large	33.0mm	16mm	0mm
Y-Large Extended reach	Large	33.0mm	26mm	6mm
Y-Large X-Extended reach	Large	33.0mm	36mm	12mm



Steerable Curve Electrophysiology Catheter Types

A - J Curve Variations


Curve Shape	Curved Length L (cm)	Curve Radius
A	3.8	1.2
B	5.1	1.6
D	6.4	2.0
F	7.6	2.4
G	8.9	2.8
J	10.2	3.3



Steerable Curve Electrophysiology Catheter Types


P - R Curve Variations

Curve Shape	Curved Length L (cm)	Curve Radius
P	12.7	2.7
Q	15.2	3.2
R	17.8	3.8




Steerable Loop Curve Electrophysiology Catheter Types

Available in the following configurations:




Curve Shape	Curved Length L (cm)	Curve Radius (cm)
B	5.1	1.6
D	6.4	2.0



Khelix Fixed Curve Electrophysiology Catheters

Quadpolar fixed curve catheters



Compatible with
X-D10 Extension
Cable



Quadpolar 2-5-2 Electrode Spacing

Ordering Information

Product Reference Number	Usable Length	Electrode Spacing	Shape Curve	Curve Size	Curve Reach	Shaft Size
5SSF-ANN-4A252	115	2-5-2	A	Nominal	Nominal	6F
5SSF-ALN-4A252	115	2-5-2	A	Large	Nominal	6F
5SSF-ANE-4A252	115	2-5-2	A	Nominal	Extended	6F
5SSF-ALE-4A252	115	2-5-2	A	Large	Extended	6F
5SSF-DNN-4A252	115	2-5-2	D	Nominal	Nominal	6F
5SSF-DLN-4A252	115	2-5-2	D	Large	Nominal	6F
5SSF-DNE-4A252	115	2-5-2	D	Nominal	Extended	6F
5SSF-DLE-4A252	115	2-5-2	D	Large	Extended	6F
5SSF-FNN-4A252	115	2-5-2	F	Nominal	Nominal	6F
5SSF-FLN-4A252	115	2-5-2	F	Large	Nominal	6F
5SSF-FNE-4A252	115	2-5-2	F	Nominal	Extended	6F
5SSF-FLE-4A252	115	2-5-2	F	Large	Extended	6F
5SSF-YNN-4A252	115	2-5-2	Y	Nominal	Nominal	6F
5SSF-YLN-4A252	115	2-5-2	Y	Large	Nominal	6F
5SSF-YNE-4A252	115	2-5-2	Y	Nominal	Extended	6F
5SSF-YLE-4A252	115	2-5-2	Y	Large	Extended	6F
5SSF-YNXE-4A252	115	2-5-2	Y	Nominal	Extra Extended	6F
5SSF-YLXE-4A252	115	2-5-2	Y	Large	Extra Extended	6F

Khelix Fixed Curve Electrophysiology Catheters

Quadpolar fixed curve catheters



Compatible with
X-D10 Extension
Cable



Quadpolar 5-5-5 Electrode Spacing

Ordering Information

Product Reference Number	Usable Length	Electrode Spacing	Shape Curve	Curve Size	Curve Reach	Shaft Size
5SSF-ANN-4A555	115	5-5-5	A	Nominal	Nominal	6F
5SSF-ALN-4A555	115	5-5-5	A	Large	Nominal	6F
5SSF-ANE-4A555	115	5-5-5	A	Nominal	Extended	6F
5SSF-ALE-4A555	115	5-5-5	A	Large	Extended	6F
5SSF-DNN-4A555	115	5-5-5	D	Nominal	Nominal	6F
5SSF-DLN-4A555	115	5-5-5	D	Large	Nominal	6F
5SSF-DNE-4A555	115	5-5-5	D	Nominal	Extended	6F
5SSF-DLE-4A555	115	5-5-5	D	Large	Extended	6F
5SSF-FNN-4A555	115	5-5-5	F	Nominal	Nominal	6F
5SSF-FLN-4A555	115	5-5-5	F	Large	Nominal	6F
5SSF-FNE-4A555	115	5-5-5	F	Nominal	Extended	6F
5SSF-FLE-4A555	115	5-5-5	F	Large	Extended	6F
5SSF-YNN-4A555	115	5-5-5	Y	Nominal	Nominal	6F
5SSF-YLN-4A555	115	5-5-5	Y	Large	Nominal	6F
5SSF-YNE-4A555	115	5-5-5	Y	Nominal	Extended	6F
5SSF-YLE-4A555	115	5-5-5	Y	Large	Extended	6F
5SSF-YNXE-4A555	115	5-5-5	Y	Nominal	Extra Extended	6F
5SSF-YLXE-4A555	115	5-5-5	Y	Large	Extra Extended	6F

Khelix Fixed Curve Electrophysiology Catheters

Quadpolar fixed curve catheters



Compatible with
X-D10 Extension
Cable



Quadpolar 10-10-10 Electrode Spacing

Ordering Information


Product Reference Number	Usable Length	Electrode Spacing	Shape Curve	Curve Size	Curve Reach	Shaft Size
5SSF-ANN-4A101010	115	10-10-10	A	Nominal	Nominal	6F
5SSF-ALN-4A101010	115	10-10-10	A	Large	Nominal	6F
5SSF-ANE-4A101010	115	10-10-10	A	Nominal	Extended	6F
5SSF-ALE-4A101010	115	10-10-10	A	Large	Extended	6F
5SSF-DNN-4A101010	115	10-10-10	D	Nominal	Nominal	6F
5SSF-DLN-4A101010	115	10-10-10	D	Large	Nominal	6F
5SSF-DNE-4A101010	115	10-10-10	D	Nominal	Extended	6F
5SSF-DLE-4A101010	115	10-10-10	D	Large	Extended	6F
5SSF-FNN-4A101010	115	10-10-10	F	Nominal	Nominal	6F
5SSF-FLN-4A101010	115	10-10-10	F	Large	Nominal	6F
5SSF-FNE-4A101010	115	10-10-10	F	Nominal	Extended	6F
5SSF-FLE-4A101010	115	10-10-10	F	Large	Extended	6F
5SSF-YNN-4A101010	115	10-10-10	Y	Nominal	Nominal	6F
5SSF-YLN-4A101010	115	10-10-10	Y	Large	Nominal	6F
5SSF-YNE-4A101010	115	10-10-10	Y	Nominal	Extended	6F
5SSF-YLE-4A101010	115	10-10-10	Y	Large	Extended	6F
5SSF-YNXE-4A101010	115	10-10-10	Y	Nominal	Extra Extended	6F
5SSF-YLXE-4A101010	115	10-10-10	Y	Large	Extra Extended	6F

Khelix Fixed Curve Electrophysiology Catheters

Hexpolar fixed curve catheters



Compatible with
X-D10 Extension
Cable



Ordering Information

Product Reference Number	Usable Length	Electrode Spacing	Shape Curve	Curve Size	Curve Reach	Shaft Size
5SSF-ANN-6A252	115	2-5-2	A	Nominal	Nominal	6F
5SSF-ALN-6A252	115	2-5-2	A	Large	Nominal	6F
5SSF-ANE-6A252	115	2-5-2	A	Nominal	Extended	6F
5SSF-ALE-6A252	115	2-5-2	A	Large	Extended	6F
5SSF-DNN-6A252	115	2-5-2	D	Nominal	Nominal	6F
5SSF-DLN-6A252	115	2-5-2	D	Large	Nominal	6F
5SSF-DNE-6A252	115	2-5-2	D	Nominal	Extended	6F
5SSF-DLE-6A252	115	2-5-2	D	Large	Extended	6F
5SSF-FNN-6A252	115	2-5-2	F	Nominal	Nominal	6F
5SSF-FLN-6A252	115	2-5-2	F	Large	Nominal	6F
5SSF-FNE-6A252	115	2-5-2	F	Nominal	Extended	6F
5SSF-FLE-6A252	115	2-5-2	F	Large	Extended	6F
5SSF-YNN-6A252	115	2-5-2	Y	Nominal	Nominal	6F
5SSF-YLN-6A252	115	2-5-2	Y	Large	Nominal	6F
5SSF-YNE-6A252	115	2-5-2	Y	Nominal	Extended	6F
5SSF-YLE-6A252	115	2-5-2	Y	Large	Extended	6F
5SSF-YNXE-6A252	115	2-5-2	Y	Nominal	Extra Extended	6F
5SSF-YLXE-6A252	115	2-5-2	Y	Large	Extra Extended	6F

Khelix Fixed Curve Electrophysiology Catheters

Hexpolar fixed curve catheters



Compatible with
X-D10 Extension
Cable



Hexpolar 5-5-5 Electrode Spacing

Ordering Information

Product Reference Number	Usable Length	Electrode Spacing	Shape Curve	Curve Size	Curve Reach	
5SSF-ANN-6A555	115	5-5-5	A	Nominal	Nominal	6F
5SSF-ALN-6A555	115	5-5-5	A	Large	Nominal	6F
5SSF-ANE-6A555	115	5-5-5	A	Nominal	Extended	6F
5SSF-ALE-6A555	115	5-5-5	A	Large	Extended	6F
5SSF-DNN-6A555	115	5-5-5	D	Nominal	Nominal	6F
5SSF-DLN-6A555	115	5-5-5	D	Large	Nominal	6F
5SSF-DNE-6A555	115	5-5-5	D	Nominal	Extended	6F
5SSF-DLE-6A555	115	5-5-5	D	Large	Extended	6F
5SSF-FNN-6A555	115	5-5-5	F	Nominal	Nominal	6F
5SSF-FLN-6A555	115	5-5-5	F	Large	Nominal	6F
5SSF-FNE-6A555	115	5-5-5	F	Nominal	Extended	6F
5SSF-FLE-6A555	115	5-5-5	F	Large	Extended	6F
5SSF-YNN-6A555	115	5-5-5	Y	Nominal	Nominal	6F
5SSF-YLN-6A555	115	5-5-5	Y	Large	Nominal	6F
5SSF-YNE-6A555	115	5-5-5	Y	Nominal	Extended	6F
5SSF-YLE-6A555	115	5-5-5	Y	Large	Extended	6F
5SSF-YNXE-6A555	115	5-5-5	Y	Nominal	Extra Extended	6F
5SSF-YLXE-6A555	115	5-5-5	Y	Large	Extra Extended	6F

Khelix Fixed Curve Electrophysiology Catheters

Hexpolar fixed curve catheters



Compatible with
X-D10 Extension
Cable



Hexpolar 10-10-10 Electrode Spacing

Ordering Information


Product Reference Number	Usable Length	Electrode Spacing	Shape Curve	Curve Size	Curve Reach	
5SSF-ANN-6A101010	115	10-10-10	A	Nominal	Nominal	6F
5SSF-ALN-6A101010	115	10-10-10	A	Large	Nominal	6F
5SSF-ANE-6A101010	115	10-10-10	A	Nominal	Extended	6F
5SSF-ALE-6A101010	115	10-10-10	A	Large	Extended	6F
5SSF-DNN-6A101010	115	10-10-10	D	Nominal	Nominal	6F
5SSF-DLN-6A101010	115	10-10-10	D	Large	Nominal	6F
5SSF-DNE-6A101010	115	10-10-10	D	Nominal	Extended	6F
5SSF-DLE-6A101010	115	10-10-10	D	Large	Extended	6F
5SSF-FNN-6A101010	115	10-10-10	F	Nominal	Nominal	6F
5SSF-FLN-6A101010	115	10-10-10	F	Large	Nominal	6F
5SSF-FNE-6A101010	115	10-10-10	F	Nominal	Extended	6F
5SSF-FLE-6A101010	115	10-10-10	F	Large	Extended	6F
5SSF-YNN-6A101010	115	10-10-10	Y	Nominal	Nominal	6F
5SSF-YLN-6A101010	115	10-10-10	Y	Large	Nominal	6F
5SSF-YNE-6A101010	115	10-10-10	Y	Nominal	Extended	6F
5SSF-YLE-6A101010	115	10-10-10	Y	Large	Extended	6F
5SSF-YNXE-6A101010	115	10-10-10	Y	Nominal	Extra Extended	6F
5SSF-YLXE-6A101010	115	10-10-10	Y	Large	Extra Extended	6F

Steerable Electrophysiology Catheters

Active Tip D curve Catheters



Compatible with
X-D10 Extension
Cable



Ordering Information

Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
SS1 - D0U - 04A6222	4	2-2-2	Active	115	6F
SS1 - D0U - 04A6252	4	2-5-2	Active	115	6F
SS1 - D0U - 04A6282	4	2-8-2	Active	115	6F
SS1 - D0U - 04A6555	4	5-5-5	Active	115	6F
SS1 - D0U - 04A62102	4	2-10-2	Active	115	6F
SS1 - D0U - 06A6222	6	2-2-2	Active	115	6F
SS1 - D0U - 06A6252	6	2-5-2	Active	115	6F
SS1 - D0U - 06A6282	6	2-8-2	Active	115	6F
SS1 - D0U - 06A6555	6	5-5-5	Active	115	6F
SS1 - D0U - 06A62102	6	2-10-2	Active	115	6F
SS1 - D0U - 08A6222	8	2-2-2	Active	115	6F
SS1 - D0U - 08A6252	8	2-5-2	Active	115	6F
SS1 - D0U - 08A6282	8	2-8-2	Active	115	6F
SS1 - D0U - 08A6555	8	5-5-5	Active	115	6F
SS1 - D0U - 08A62102	8	2-10-2	Active	115	6F
SS1 - D0U - 10A6222	10	2-2-2	Active	115	6F
SS1 - D0U - 10A6252	10	2-5-2	Active	115	6F
SS1 - D0U - 10A6282	10	2-8-2	Active	115	6F
SS1 - D0U - 10A6555	10	5-5-5	Active	115	6F
SS1 - D0U - 10A62102	10	2-10-2	Active	115	6F

Steerable Electrophysiology Catheters

Passive Tip D curve Catheters



Compatible with
X-D10 Extension
Cable




Ordering Information


Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
SS1 - D0U - 04P6222	4	2-2-2	Passive	115	6F
SS1 - D0U - 04P6252	4	2-5-2	Passive	115	6F
SS1 - D0U - 04P6282	4	2-8-2	Passive	115	6F
SS1 - D0U - 04P6555	4	5-5-5	Passive	115	6F
SS1 - D0U - 04P62102	4	2-10-2	Passive	115	6F
SS1 - D0U - 06P6222	6	2-2-2	Passive	115	6F
SS1 - D0U - 06P6252	6	2-5-2	Passive	115	6F
SS1 - D0U - 06P6282	6	2-8-2	Passive	115	6F
SS1 - D0U - 06P6555	6	5-5-5	Passive	115	6F
SS1 - D0U - 06P62102	6	2-10-2	Passive	115	6F
SS1 - D0U - 08P6222	8	2-2-2	Passive	115	6F
SS1 - D0U - 08P6252	8	2-5-2	Passive	115	6F
SS1 - D0U - 08P6282	8	2-8-2	Passive	115	6F
SS1 - D0U - 08P6555	8	5-5-5	Passive	115	6F
SS1 - D0U - 08P62102	8	2-10-2	Passive	115	6F
SS1 - D0U - 10P6222	10	2-2-2	Passive	115	6F
SS1 - D0U - 10P6252	10	2-5-2	Passive	115	6F
SS1 - D0U - 10P6282	10	2-8-2	Passive	115	6F
SS1 - D0U - 10P6555	10	5-5-5	Passive	115	6F
SS1 - D0U - 10P62102	10	2-10-2	Passive	115	6F

Steerable Electrophysiology Catheters

Active Tip F curve Catheters



Compatible with
X-D10 Extension
Cable




Ordering Information

Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
SS1 - F0U - 04A6222	4	2-2-2	Active	115	6F
SS1 - F0U - 04A6252	4	2-5-2	Active	115	6F
SS1 - F0U - 04A6282	4	2-8-2	Active	115	6F
SS1 - F0U - 04A6555	4	5-5-5	Active	115	6F
SS1 - F0U - 04A62102	4	2-10-2	Active	115	6F
SS1 - F0U - 06A6222	6	2-2-2	Active	115	6F
SS1 - F0U - 06A6252	6	2-5-2	Active	115	6F
SS1 - F0U - 06A6282	6	2-8-2	Active	115	6F
SS1 - F0U - 06A6555	6	5-5-5	Active	115	6F
SS1 - F0U - 06A62102	6	2-10-2	Active	115	6F
SS1 - F0U - 08A6222	8	2-2-2	Active	115	6F
SS1 - F0U - 08A6252	8	2-5-2	Active	115	6F
SS1 - F0U - 08A6282	8	2-8-2	Active	115	6F
SS1 - F0U - 08A6555	8	5-5-5	Active	115	6F
SS1 - F0U - 08A62102	8	2-10-2	Active	115	6F
SS1 - F0U - 10A6222	10	2-2-2	Active	115	6F
SS1 - F0U - 10A6252	10	2-5-2	Active	115	6F
SS1 - F0U - 10A6282	10	2-8-2	Active	115	6F
SS1 - F0U - 10A6555	10	5-5-5	Active	115	6F
SS1 - F0U - 10A62102	10	2-10-2	Active	115	6F

Steerable Electrophysiology Catheters

Passive Tip F curve Catheters



Compatible with
X-D10 Extension
Cable




Ordering Information


Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
SS1 - F0U - 04P6222	4	2-2-2	Passive	115	6F
SS1 - F0U - 04P6252	4	2-5-2	Passive	115	6F
SS1 - F0U - 04P6282	4	2-8-2	Passive	115	6F
SS1 - F0U - 04P6555	4	5-5-5	Passive	115	6F
SS1 - F0U - 04P62102	4	2-10-2	Passive	115	6F
SS1 - F0U - 06P6222	6	2-2-2	Passive	115	6F
SS1 - F0U - 06P6252	6	2-5-2	Passive	115	6F
SS1 - F0U - 06P6282	6	2-8-2	Passive	115	6F
SS1 - F0U - 06P6555	6	5-5-5	Passive	115	6F
SS1 - F0U - 06P62102	6	2-10-2	Passive	115	6F
SS1 - F0U - 08P6222	8	2-2-2	Passive	115	6F
SS1 - F0U - 08P6252	8	2-5-2	Passive	115	6F
SS1 - F0U - 08P6282	8	2-8-2	Passive	115	6F
SS1 - F0U - 08P6555	8	5-5-5	Passive	115	6F
SS1 - F0U - 08P62102	8	2-10-2	Passive	115	6F
SS1 - F0U - 10P6222	10	2-2-2	Passive	115	6F
SS1 - F0U - 10P6252	10	2-5-2	Passive	115	6F
SS1 - F0U - 10P6282	10	2-8-2	Passive	115	6F
SS1 - F0U - 10P6555	10	5-5-5	Passive	115	6F
SS1 - F0U - 10P62102	10	2-10-2	Passive	115	6F

Steerable Electrophysiology Catheters

Active Tip G curve Catheters



Compatible with
X-D10 Extension
Cable



Ordering Information


Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
SS1 - G0U - 04A6222	4	2-2-2	Active	115	6F
SS1 - G0U - 04A6252	4	2-5-2	Active	115	6F
SS1 - G0U - 04A6282	4	2-8-2	Active	115	6F
SS1 - G0U - 04A6555	4	5-5-5	Active	115	6F
SS1 - G0U - 04A62102	4	2-10-2	Active	115	6F
SS1 - G0U - 06A6222	6	2-2-2	Active	115	6F
SS1 - G0U - 06A6252	6	2-5-2	Active	115	6F
SS1 - G0U - 06A6282	6	2-8-2	Active	115	6F
SS1 - G0U - 06A6555	6	5-5-5	Active	115	6F
SS1 - G0U - 06A62102	6	2-10-2	Active	115	6F
SS1 - G0U - 08A6222	8	2-2-2	Active	115	6F
SS1 - G0U - 08A6252	8	2-5-2	Active	115	6F
SS1 - G0U - 08A6282	8	2-8-2	Active	115	6F
SS1 - G0U - 08A6555	8	5-5-5	Active	115	6F
SS1 - G0U - 08A62102	8	2-10-2	Active	115	6F
SS1 - G0U - 10A6222	10	2-2-2	Active	115	6F
SS1 - G0U - 10A6252	10	2-5-2	Active	115	6F
SS1 - G0U - 10A6282	10	2-8-2	Active	115	6F
SS1 - G0U - 10A6555	10	5-5-5	Active	115	6F
SS1 - G0U - 10A62102	10	2-10-2	Active	115	6F

Steerable Electrophysiology Catheters

Passive Tip G curve Catheters



Compatible with
X-D10 Extension
Cable




Ordering Information

Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
SS1 - G0U - 04P6222	4	2-2-2	Passive	115	6F
SS1 - G0U - 04P6252	4	2-5-2	Passive	115	6F
SS1 - G0U - 04P6282	4	2-8-2	Passive	115	6F
SS1 - G0U - 04P6555	4	5-5-5	Passive	115	6F
SS1 - G0U - 04P62102	4	2-10-2	Passive	115	6F
SS1 - G0U - 06P6222	6	2-2-2	Passive	115	6F
SS1 - G0U - 06P6252	6	2-5-2	Passive	115	6F
SS1 - G0U - 06P6282	6	2-8-2	Passive	115	6F
SS1 - G0U - 06P6555	6	5-5-5	Passive	115	6F
SS1 - G0U - 06P62102	6	2-10-2	Passive	115	6F
SS1 - G0U - 08P6222	8	2-2-2	Passive	115	6F
SS1 - G0U - 08P6252	8	2-5-2	Passive	115	6F
SS1 - G0U - 08P6282	8	2-8-2	Passive	115	6F
SS1 - G0U - 08P6555	8	5-5-5	Passive	115	6F
SS1 - G0U - 08P62102	8	2-10-2	Passive	115	6F
SS1 - G0U - 10P6222	10	2-2-2	Passive	115	6F
SS1 - G0U - 10P6252	10	2-5-2	Passive	115	6F
SS1 - G0U - 10P6282	10	2-8-2	Passive	115	6F
SS1 - G0U - 10P6555	10	5-5-5	Passive	115	6F
SS1 - G0U - 10P62102	10	2-10-2	Passive	115	6F

Steerable Electrophysiology Catheters

Active Tip J curve Catheters



Compatible with
X-D10 Extension
Cable




Ordering Information

Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
SS1 - J0U - 04A6222	4	2-2-2	Active	115	6F
SS1 - J0U - 04A6252	4	2-5-2	Active	115	6F
SS1 - J0U - 04A6282	4	2-8-2	Active	115	6F
SS1 - J0U - 04A6555	4	5-5-5	Active	115	6F
SS1 - J0U - 04A62102	4	2-10-2	Active	115	6F
SS1 - J0U - 06A6222	6	2-2-2	Active	115	6F
SS1 - J0U - 06A6252	6	2-5-2	Active	115	6F
SS1 - J0U - 06A6282	6	2-8-2	Active	115	6F
SS1 - J0U - 06A6555	6	5-5-5	Active	115	6F
SS1 - J0U - 06A62102	6	2-10-2	Active	115	6F
SS1 - J0U - 08A6222	8	2-2-2	Active	115	6F
SS1 - J0U - 08A6252	8	2-5-2	Active	115	6F
SS1 - J0U - 08A6282	8	2-8-2	Active	115	6F
SS1 - J0U - 08A6555	8	5-5-5	Active	115	6F
SS1 - J0U - 08A62102	8	2-10-2	Active	115	6F
SS1 - J0U - 10A6222	10	2-2-2	Active	115	6F
SS1 - J0U - 10A6252	10	2-5-2	Active	115	6F
SS1 - J0U - 10A6282	10	2-8-2	Active	115	6F
SS1 - J0U - 10A6555	10	5-5-5	Active	115	6F
SS1 - J0U - 10A62102	10	2-10-2	Active	115	6F

Steerable Electrophysiology Catheters

Passive Tip J curve Catheters



Compatible with
X-D10 Extension
Cable




Ordering Information


Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
SS1 - J0U - 04P6222	4	2-2-2	Passive	115	6F
SS1 - J0U - 04P6252	4	2-5-2	Passive	115	6F
SS1 - J0U - 04P6282	4	2-8-2	Passive	115	6F
SS1 - J0U - 04P6555	4	5-5-5	Passive	115	6F
SS1 - J0U - 04P62102	4	2-10-2	Passive	115	6F
SS1 - J0U - 06P6222	6	2-2-2	Passive	115	6F
SS1 - J0U - 06P6252	6	2-5-2	Passive	115	6F
SS1 - J0U - 06P6282	6	2-8-2	Passive	115	6F
SS1 - J0U - 06P6555	6	5-5-5	Passive	115	6F
SS1 - J0U - 06P62102	6	2-10-2	Passive	115	6F
SS1 - J0U - 08P6222	8	2-2-2	Passive	115	6F
SS1 - J0U - 08P6252	8	2-5-2	Passive	115	6F
SS1 - J0U - 08P6282	8	2-8-2	Passive	115	6F
SS1 - J0U - 08P6555	8	5-5-5	Passive	115	6F
SS1 - J0U - 08P62102	8	2-10-2	Passive	115	6F
SS1 - J0U - 10P6222	10	2-2-2	Passive	115	6F
SS1 - J0U - 10P6252	10	2-5-2	Passive	115	6F
SS1 - J0U - 10P6282	10	2-8-2	Passive	115	6F
SS1 - J0U - 10P6555	10	5-5-5	Passive	115	6F
SS1 - J0U - 10P62102	10	2-10-2	Passive	115	6F

Steerable Electrophysiology Catheters

Active Tip P curve Catheters



Compatible with
X-D10 Extension
Cable



Ordering Information


Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
SS1 - P0U - 04A6222	4	2-2-2	Active	115	6F
SS1 - P0U - 04A6252	4	2-5-2	Active	115	6F
SS1 - P0U - 04A6282	4	2-8-2	Active	115	6F
SS1 - P0U - 04A6555	4	5-5-5	Active	115	6F
SS1 - P0U - 04A62102	4	2-10-2	Active	115	6F
SS1 - P0U - 04A6101010	4	10-10-10	Active	115	6F
SS1 - P0U - 06A6222	6	2-2-2	Active	115	6F
SS1 - P0U - 06A6252	6	2-5-2	Active	115	6F
SS1 - P0U - 06A6282	6	2-8-2	Active	115	6F
SS1 - P0U - 06A6555	6	5-5-5	Active	115	6F
SS1 - P0U - 06A62102	6	2-10-2	Active	115	6F
SS1 - P0U - 06A6101010	6	10-10-10	Active	115	6F
SS1 - P0U - 08A6222	8	2-2-2	Active	115	6F
SS1 - P0U - 08A6252	8	2-5-2	Active	115	6F
SS1 - P0U - 08A6282	8	2-8-2	Active	115	6F
SS1 - P0U - 08A6555	8	5-5-5	Active	115	6F
SS1 - P0U - 08A62102	8	2-10-2	Active	115	6F
SS1 - P0U - 08A6101010	8	10-10-10	Active	115	6F
SS1 - P0U - 10A6222	10	2-2-2	Active	115	6F
SS1 - P0U - 10A6252	10	2-5-2	Active	115	6F
SS1 - P0U - 10A6282	10	2-8-2	Active	115	6F
SS1 - P0U - 10A6555	10	5-5-5	Active	115	6F
SS1 - P0U - 10A62102	10	2-10-2	Active	115	6F
SS1 - P0U - 10A6101010	10	10-10-10	Active	115	6F

Steerable Electrophysiology Catheters

Passive Tip P curve Catheters



Compatible with
X-D10 Extension
Cable



Ordering Information


Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
SS1 - P0U - 04P6222	4	2-2-2	Passive	115	6F
SS1 - P0U - 04P6252	4	2-5-2	Passive	115	6F
SS1 - P0U - 04P6282	4	2-8-2	Passive	115	6F
SS1 - P0U - 04P6555	4	5-5-5	Passive	115	6F
SS1 - P0U - 04P62102	4	2-10-2	Passive	115	6F
SS1 - P0U - 04P6101010	4	10-10-10	Passive	115	6F
SS1 - P0U - 06P6222	6	2-2-2	Passive	115	6F
SS1 - P0U - 06P6252	6	2-5-2	Passive	115	6F
SS1 - P0U - 06P6282	6	2-8-2	Passive	115	6F
SS1 - P0U - 06P6555	6	5-5-5	Passive	115	6F
SS1 - P0U - 06P62102	6	2-10-2	Passive	115	6F
SS1 - P0U - 06P6101010	6	10-10-10	Passive	115	6F
SS1 - P0U - 08P6222	8	2-2-2	Passive	115	6F
SS1 - P0U - 08P6252	8	2-5-2	Passive	115	6F
SS1 - P0U - 08P6282	8	2-8-2	Passive	115	6F
SS1 - P0U - 08P6555	8	5-5-5	Passive	115	6F
SS1 - P0U - 08P62102	8	2-10-2	Passive	115	6F
SS1 - P0U - 08P6101010	8	10-10-10	Passive	115	6F
SS1 - P0U - 10P6222	10	2-2-2	Passive	115	6F
SS1 - P0U - 10P6252	10	2-5-2	Passive	115	6F
SS1 - P0U - 10P6282	10	2-8-2	Passive	115	6F
SS1 - P0U - 10P6555	10	5-5-5	Passive	115	6F
SS1 - P0U - 10P62102	10	2-10-2	Passive	115	6F
SS1 - P0U - 10P6101010	10	10-10-10	Passive	115	6F

Steerable Electrophysiology Catheters

Active Tip Q curve Catheters



Compatible with
X-D10 Extension
Cable



Ordering Information


Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
SS1 - Q0U - 04A6222	4	2-2-2	Active	115	6F
SS1 - Q0U - 04A6252	4	2-5-2	Active	115	6F
SS1 - Q0U - 04A6282	4	2-8-2	Active	115	6F
SS1 - Q0U - 04A6555	4	5-5-5	Active	115	6F
SS1 - Q0U - 04A62102	4	2-10-2	Active	115	6F
SS1 - Q0U - 04A6101010	4	10-10-10	Active	115	6F
SS1 - Q0U - 06A6222	6	2-2-2	Active	115	6F
SS1 - Q0U - 06A6252	6	2-5-2	Active	115	6F
SS1 - Q0U - 06A6282	6	2-8-2	Active	115	6F
SS1 - Q0U - 06A6555	6	5-5-5	Active	115	6F
SS1 - Q0U - 06A62102	6	2-10-2	Active	115	6F
SS1 - Q0U - 06A6101010	6	10-10-10	Active	115	6F
SS1 - Q0U - 08A6222	8	2-2-2	Active	115	6F
SS1 - Q0U - 08A6252	8	2-5-2	Active	115	6F
SS1 - Q0U - 08A6282	8	2-8-2	Active	115	6F
SS1 - Q0U - 08A6555	8	5-5-5	Active	115	6F
SS1 - Q0U - 08A62102	8	2-10-2	Active	115	6F
SS1 - Q0U - 08A6101010	8	10-10-10	Active	115	6F
SS1 - Q0U - 10A6222	10	2-2-2	Active	115	6F
SS1 - Q0U - 10A6252	10	2-5-2	Active	115	6F
SS1 - Q0U - 10A6282	10	2-8-2	Active	115	6F
SS1 - Q0U - 10A6555	10	5-5-5	Active	115	6F
SS1 - Q0U - 10A62102	10	2-10-2	Active	115	6F
SS1 - Q0U - 10A6101010	10	10-10-10	Active	115	6F

Steerable Electrophysiology Catheters

Passive Tip Q curve Catheters



Compatible with
X-D10 Extension
Cable



Ordering Information


Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
SS1 - Q0U - 04P6222	4	2-2-2	Passive	115	6F
SS1 - Q0U - 04P6252	4	2-5-2	Passive	115	6F
SS1 - Q0U - 04P6282	4	2-8-2	Passive	115	6F
SS1 - Q0U - 04P6555	4	5-5-5	Passive	115	6F
SS1 - Q0U - 04P62102	4	2-10-2	Passive	115	6F
SS1 - Q0U - 04P6101010	4	10-10-10	Passive	115	6F
SS1 - Q0U - 06P6222	6	2-2-2	Passive	115	6F
SS1 - Q0U - 06P6252	6	2-5-2	Passive	115	6F
SS1 - Q0U - 06P6282	6	2-8-2	Passive	115	6F
SS1 - Q0U - 06P6555	6	5-5-5	Passive	115	6F
SS1 - Q0U - 06P62102	6	2-10-2	Passive	115	6F
SS1 - Q0U - 06P6101010	6	10-10-10	Passive	115	6F
SS1 - Q0U - 08P6222	8	2-2-2	Passive	115	6F
SS1 - Q0U - 08P6252	8	2-5-2	Passive	115	6F
SS1 - Q0U - 08P6282	8	2-8-2	Passive	115	6F
SS1 - Q0U - 08P6555	8	5-5-5	Passive	115	6F
SS1 - Q0U - 08P62102	8	2-10-2	Passive	115	6F
SS1 - Q0U - 08P6101010	8	10-10-10	Passive	115	6F
SS1 - Q0U - 10P6222	10	2-2-2	Passive	115	6F
SS1 - Q0U - 10P6252	10	2-5-2	Passive	115	6F
SS1 - Q0U - 10P6282	10	2-8-2	Passive	115	6F
SS1 - Q0U - 10P6555	10	5-5-5	Passive	115	6F
SS1 - Q0U - 10P62102	10	2-10-2	Passive	115	6F
SS1 - Q0U - 10P6101010	10	10-10-10	Passive	115	6F

Steerable Electrophysiology Catheters

Active Tip R curve Catheters



Compatible with
X-D10 Extension
Cable




Ordering Information


Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
SS1 - R0U - 04A6222	4	2-2-2	Active	115	6F
SS1 - R0U - 04A6252	4	2-5-2	Active	115	6F
SS1 - R0U - 04A6282	4	2-8-2	Active	115	6F
SS1 - R0U - 04A6555	4	5-5-5	Active	115	6F
SS1 - R0U - 04A62102	4	2-10-2	Active	115	6F
SS1 - R0U - 04A6101010	4	10-10-10	Active	115	6F
SS1 - R0U - 06A6222	6	2-2-2	Active	115	6F
SS1 - R0U - 06A6252	6	2-5-2	Active	115	6F
SS1 - R0U - 06A6282	6	2-8-2	Active	115	6F
SS1 - R0U - 06A6555	6	5-5-5	Active	115	6F
SS1 - R0U - 06A62102	6	2-10-2	Active	115	6F
SS1 - R0U - 06A6101010	6	10-10-10	Active	115	6F
SS1 - R0U - 08A6222	8	2-2-2	Active	115	6F
SS1 - R0U - 08A6252	8	2-5-2	Active	115	6F
SS1 - R0U - 08A6282	8	2-8-2	Active	115	6F
SS1 - R0U - 08A6555	8	5-5-5	Active	115	6F
SS1 - R0U - 08A62102	8	2-10-2	Active	115	6F
SS1 - R0U - 08A6101010	8	10-10-10	Active	115	6F
SS1 - R0U - 10A6222	10	2-2-2	Active	115	6F
SS1 - R0U - 10A6252	10	2-5-2	Active	115	6F
SS1 - R0U - 10A6282	10	2-8-2	Active	115	6F
SS1 - R0U - 10A6555	10	5-5-5	Active	115	6F
SS1 - R0U - 10A62102	10	2-10-2	Active	115	6F
SS1 - R0U - 10A6101010	10	10-10-10	Active	115	6F

Steerable Electrophysiology Catheters

Passive Tip R curve Catheters



Compatible with
X-D10 Extension
Cable



Ordering Information


Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
SS1 - R0U - 04P6222	4	2-2-2	Passive	115	6F
SS1 - R0U - 04P6252	4	2-5-2	Passive	115	6F
SS1 - R0U - 04P6282	4	2-8-2	Passive	115	6F
SS1 - R0U - 04P6555	4	5-5-5	Passive	115	6F
SS1 - R0U - 04P62102	4	2-10-2	Passive	115	6F
SS1 - R0U - 04P6101010	4	10-10-10	Passive	115	6F
SS1 - R0U - 06P6222	6	2-2-2	Passive	115	6F
SS1 - R0U - 06P6252	6	2-5-2	Passive	115	6F
SS1 - R0U - 06P6282	6	2-8-2	Passive	115	6F
SS1 - R0U - 06P6555	6	5-5-5	Passive	115	6F
SS1 - R0U - 06P62102	6	2-10-2	Passive	115	6F
SS1 - R0U - 06P6101010	6	10-10-10	Passive	115	6F
SS1 - R0U - 08P6222	8	2-2-2	Passive	115	6F
SS1 - R0U - 08P6252	8	2-5-2	Passive	115	6F
SS1 - R0U - 08P6282	8	2-8-2	Passive	115	6F
SS1 - R0U - 08P6555	8	5-5-5	Passive	115	6F
SS1 - R0U - 08P62102	8	2-10-2	Passive	115	6F
SS1 - R0U - 08P6101010	8	10-10-10	Passive	115	6F
SS1 - R0U - 10P6222	10	2-2-2	Passive	115	6F
SS1 - R0U - 10P6252	10	2-5-2	Passive	115	6F
SS1 - R0U - 10P6282	10	2-8-2	Passive	115	6F
SS1 - R0U - 10P6555	10	5-5-5	Passive	115	6F
SS1 - R0U - 10P62102	10	2-10-2	Passive	115	6F
SS1 - R0U - 10P6101010	10	10-10-10	Passive	115	6F

Steerable Electrophysiology Catheters

Active Tip Q curve Catheters



Compatible with
X-D10 Extension
Cable



Ordering Information


Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
SS1 - Q2U - 04A6222	4	2-2-2	Active	115	6F
SS1 - Q2U - 04A6252	4	2-5-2	Active	115	6F
SS1 - Q2U - 04A6282	4	2-8-2	Active	115	6F
SS1 - Q2U - 04A6555	4	5-5-5	Active	115	6F
SS1 - Q2U - 04A62102	4	2-10-2	Active	115	6F
SS1 - Q2U - 04A6101010	4	10-10-10	Active	115	6F
SS1 - Q2U - 06A6222	6	2-2-2	Active	115	6F
SS1 - Q2U - 06A6252	6	2-5-2	Active	115	6F
SS1 - Q2U - 06A6282	6	2-8-2	Active	115	6F
SS1 - Q2U - 06A6555	6	5-5-5	Active	115	6F
SS1 - Q2U - 06A62102	6	2-10-2	Active	115	6F
SS1 - Q2U - 06A6101010	6	10-10-10	Active	115	6F
SS1 - Q2U - 08A6222	8	2-2-2	Active	115	6F
SS1 - Q2U - 08A6252	8	2-5-2	Active	115	6F
SS1 - Q2U - 08A6282	8	2-8-2	Active	115	6F
SS1 - Q2U - 08A6555	8	5-5-5	Active	115	6F
SS1 - Q2U - 08A62102	8	2-10-2	Active	115	6F
SS1 - Q2U - 08A6101010	8	10-10-10	Active	115	6F
SS1 - Q2U - 10A6222	10	2-2-2	Active	115	6F
SS1 - Q2U - 10A6252	10	2-5-2	Active	115	6F
SS1 - Q2U - 10A6282	10	2-8-2	Active	115	6F
SS1 - Q2U - 10A6555	10	5-5-5	Active	115	6F
SS1 - Q2U - 10A62102	10	2-10-2	Active	115	6F
SS1 - Q2U - 10A6101010	10	10-10-10	Active	115	6F

Steerable Electrophysiology Catheters

Passive Tip Q curve Catheters



Compatible with
X-D10 Extension
Cable



Ordering Information


Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
SS1 - Q2U - 04P6222	4	2-2-2	Passive	115	6F
SS1 - Q2U - 04P6252	4	2-5-2	Passive	115	6F
SS1 - Q2U - 04P6282	4	2-8-2	Passive	115	6F
SS1 - Q2U - 04P6555	4	5-5-5	Passive	115	6F
SS1 - Q2U - 04P62102	4	2-10-2	Passive	115	6F
SS1 - Q2U - 04P6101010	4	10-10-10	Passive	115	6F
SS1 - Q2U - 06P6222	6	2-2-2	Passive	115	6F
SS1 - Q2U - 06P6252	6	2-5-2	Passive	115	6F
SS1 - Q2U - 06P6282	6	2-8-2	Passive	115	6F
SS1 - Q2U - 06P6555	6	5-5-5	Passive	115	6F
SS1 - Q2U - 06P62102	6	2-10-2	Passive	115	6F
SS1 - Q2U - 06P6101010	6	10-10-10	Passive	115	6F
SS1 - Q2U - 08P6222	8	2-2-2	Passive	115	6F
SS1 - Q2U - 08P6252	8	2-5-2	Passive	115	6F
SS1 - Q2U - 08P6282	8	2-8-2	Passive	115	6F
SS1 - Q2U - 08P6555	8	5-5-5	Passive	115	6F
SS1 - Q2U - 08P62102	8	2-10-2	Passive	115	6F
SS1 - Q2U - 08P6101010	8	10-10-10	Passive	115	6F
SS1 - Q2U - 10P6222	10	2-2-2	Passive	115	6F
SS1 - Q2U - 10P6252	10	2-5-2	Passive	115	6F
SS1 - Q2U - 10P6282	10	2-8-2	Passive	115	6F
SS1 - Q2U - 10P6555	10	5-5-5	Passive	115	6F
SS1 - Q2U - 10P62102	10	2-10-2	Passive	115	6F
SS1 - Q2U - 10P6101010	10	10-10-10	Passive	115	6F

Steerable Electrophysiology Catheters

Active Tip R curve Catheters



Compatible with
X-D10 Extension
Cable



Ordering Information


Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
SS1 - R2U - 04A6222	4	2-2-2	Active	115	6F
SS1 - R2U - 04A6252	4	2-5-2	Active	115	6F
SS1 - R2U - 04A6282	4	2-8-2	Active	115	6F
SS1 - R2U - 04A6555	4	5-5-5	Active	115	6F
SS1 - R2U - 04A62102	4	2-10-2	Active	115	6F
SS1 - R2U - 04A6101010	4	10-10-10	Active	115	6F
SS1 - R2U - 06A6222	6	2-2-2	Active	115	6F
SS1 - R2U - 06A6252	6	2-5-2	Active	115	6F
SS1 - R2U - 06A6282	6	2-8-2	Active	115	6F
SS1 - R2U - 06A6555	6	5-5-5	Active	115	6F
SS1 - R2U - 06A62102	6	2-10-2	Active	115	6F
SS1 - R2U - 06A6101010	6	10-10-10	Active	115	6F
SS1 - R2U - 08A6222	8	2-2-2	Active	115	6F
SS1 - R2U - 08A6252	8	2-5-2	Active	115	6F
SS1 - R2U - 08A6282	8	2-8-2	Active	115	6F
SS1 - R2U - 08A6555	8	5-5-5	Active	115	6F
SS1 - R2U - 08A62102	8	2-10-2	Active	115	6F
SS1 - R2U - 08A6101010	8	10-10-10	Active	115	6F
SS1 - R2U - 10A6222	10	2-2-2	Active	115	6F
SS1 - R2U - 10A6252	10	2-5-2	Active	115	6F
SS1 - R2U - 10A6282	10	2-8-2	Active	115	6F
SS1 - R2U - 10A6555	10	5-5-5	Active	115	6F
SS1 - R2U - 10A62102	10	2-10-2	Active	115	6F
SS1 - R2U - 10A6101010	10	10-10-10	Active	115	6F

Steerable Electrophysiology Catheters

Passive Tip R curve Catheters



Compatible with
X-D10 Extension
Cable



Ordering Information


Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
SS1 - R2U - 04P6222	4	2-2-2	Passive	115	6F
SS1 - R2U - 04P6252	4	2-5-2	Passive	115	6F
SS1 - R2U - 04P6282	4	2-8-2	Passive	115	6F
SS1 - R2U - 04P6555	4	5-5-5	Passive	115	6F
SS1 - R2U - 04P62102	4	2-10-2	Passive	115	6F
SS1 - R2U - 04P6101010	4	10-10-10	Passive	115	6F
SS1 - R2U - 06P6222	6	2-2-2	Passive	115	6F
SS1 - R2U - 06P6252	6	2-5-2	Passive	115	6F
SS1 - R2U - 06P6282	6	2-8-2	Passive	115	6F
SS1 - R2U - 06P6555	6	5-5-5	Passive	115	6F
SS1 - R2U - 06P62102	6	2-10-2	Passive	115	6F
SS1 - R2U - 06P6101010	6	10-10-10	Passive	115	6F
SS1 - R2U - 08P6222	8	2-2-2	Passive	115	6F
SS1 - R2U - 08P6252	8	2-5-2	Passive	115	6F
SS1 - R2U - 08P6282	8	2-8-2	Passive	115	6F
SS1 - R2U - 08P6555	8	5-5-5	Passive	115	6F
SS1 - R2U - 08P62102	8	2-10-2	Passive	115	6F
SS1 - R2U - 08P6101010	8	10-10-10	Passive	115	6F
SS1 - R2U - 10P6222	10	2-2-2	Passive	115	6F
SS1 - R2U - 10P6252	10	2-5-2	Passive	115	6F
SS1 - R2U - 10P6282	10	2-8-2	Passive	115	6F
SS1 - R2U - 10P6555	10	5-5-5	Passive	115	6F
SS1 - R2U - 10P62102	10	2-10-2	Passive	115	6F
SS1 - R2U - 10P6101010	10	10-10-10	Passive	115	6F

Steerable Electrophysiology Catheters

P Curve Catheters



Compatible with
X-D20 Extension
Cable



Ordering Information


Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
SS2 - P0U - 20A6242	20	2-4-2	Active	115	6F
SS2 - P0U - 20A6282	20	2-8-2	Active	115	6F
SS2 - P0U - 20A62102	20	2-10-2	Active	115	6F
SS2 - P0U - 20A62122	20	2-12-2	Active	115	6F
SS2 - P0U - 20A62132	20	2-13-2	Active	115	6F
SS2 - P0U - 20P6242	20	2-4-2	Passive	115	6F
SS2 - P0U - 20P6282	20	2-8-2	Passive	115	6F
SS2 - P0U - 20P62102	20	2-10-2	Passive	115	6F
SS2 - P0U - 20P62122	20	2-12-2	Passive	115	6F
SS2 - P0U - 20P62132	20	2-13-2	Passive	115	6F
SS2 - P0U - 20P620282	20	20-2-8-2	Passive	115	6F

Steerable Electrophysiology Catheters

Q Curve Catheters



Compatible with
X-D20 Extension
Cable




Ordering Information


Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
SS2 - Q0U - 20A6242	20	2-4-2	Active	115	6F
SS2 - Q0U - 20A6282	20	2-8-2	Active	115	6F
SS2 - Q0U - 20A62102	20	2-10-2	Active	115	6F
SS2 - Q0U - 20A62122	20	2-12-2	Active	115	6F
SS2 - Q0U - 20A62132	20	2-13-2	Active	115	6F
SS2 - Q0U - 20P6242	20	2-4-2	Passive	115	6F
SS2 - Q0U - 20P6282	20	2-8-2	Passive	115	6F
SS2 - Q0U - 20P62102	20	2-10-2	Passive	115	6F
SS2 - Q0U - 20P62122	20	2-12-2	Passive	115	6F
SS2 - Q0U - 20P62132	20	2-13-2	Passive	115	6F
SS2 - Q0U - 20P620282	20	20-2-8-2	Passive	115	6F

Steerable Electrophysiology Catheters

R Curve Catheters



Compatible with
X-D20 Extension
Cable





Ordering Information

Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
SS2 - R0U - 20A6242	20	2-4-2	Active	115	6F
SS2 - R0U - 20A6282	20	2-8-2	Active	115	6F
SS2 - R0U - 20A62102	20	2-10-2	Active	115	6F
SS2 - R0U - 20A62122	20	2-12-2	Active	115	6F
SS2 - R0U - 20A62132	20	2-13-2	Active	115	6F
SS2 - R0U - 20P6242	20	2-4-2	Passive	115	6F
SS2 - R0U - 20P6282	20	2-8-2	Passive	115	6F
SS2 - R0U - 20P62102	20	2-10-2	Passive	115	6F
SS2 - R0U - 20P62122	20	2-12-2	Passive	115	6F
SS2 - R0U - 20P62132	20	2-13-2	Passive	115	6F
SS2 - R0U - 20P620282	20	20-2-8-2	Passive	115	6F

Steerable Curve Catheters With Reuse Functionality

Post procedure,
unlock and
remove the
electrode sheath.



Handle and stylet
are sterilized
in CSSD via
standard steam
autoclave cycle.

A new EtO
sterilized
electrode sheath
is fitted prior
to next use.




Steerable Electrophysiology Catheters With Reuse Functionality

D Curve Catheters



Compatible with
X-D10 Extension
Cable



Ordering Information


Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
MS1 - D0U - 04P6222	4	2-2-2	Passive	115	6F
MS1 - D0U - 04P6252	4	2-5-2	Passive	115	6F
MS1 - D0U - 04P6555	4	5-5-5	Passive	115	6F
MS1 - D0U - 04P6282	4	2-8-2	Passive	115	6F
MS1 - D0U - 04P62102	4	2-10-2	Passive	115	6F
MS1 - D0U - 04P6101010	4	10-10-10	Passive	115	6F
MS1 - D0U - 06P6222	6	2-2-2	Passive	115	6F
MS1 - D0U - 06P6252	6	2-5-2	Passive	115	6F
MS1 - D0U - 06P6555	6	5-5-5	Passive	115	6F
MS1 - D0U - 06P6282	6	2-8-2	Passive	115	6F
MS1 - D0U - 06P62102	6	2-10-2	Passive	115	6F
MS1 - D0U - 06P6101010	6	10-10-10	Passive	115	6F
MS1 - D0U - 08P6222	8	2-2-2	Passive	115	6F
MS1 - D0U - 08P6252	8	2-5-2	Passive	115	6F
MS1 - D0U - 08P6555	8	5-5-5	Passive	115	6F
MS1 - D0U - 08P6282	8	2-8-2	Passive	115	6F
MS1 - D0U - 08P62102	8	2-10-2	Passive	115	6F
MS1 - D0U - 08P6101010	8	10-10-10	Passive	115	6F
MS1 - D0U - 10P6222	10	2-2-2	Passive	115	6F
MS1 - D0U - 10P6252	10	2-5-2	Passive	115	6F
MS1 - D0U - 10P6555	10	5-5-5	Passive	115	6F
MS1 - D0U - 10P6282	10	2-8-2	Passive	115	6F
MS1 - D0U - 10P62102	10	2-10-2	Passive	115	6F
MS1 - D0U - 10P6101010	10	10-10-10	Passive	115	6F

Steerable Electrophysiology Catheters With Reuse Functionality

F Curve Catheters



Compatible with
X-D10 Extension
Cable



Ordering Information


Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
MS1 - F0U - 04P6222	4	2-2-2	Passive	115	6F
MS1 - F0U - 04P6252	4	2-5-2	Passive	115	6F
MS1 - F0U - 04P6555	4	5-5-5	Passive	115	6F
MS1 - F0U - 04P6282	4	2-8-2	Passive	115	6F
MS1 - F0U - 04P62102	4	2-10-2	Passive	115	6F
MS1 - F0U - 04P6101010	4	10-10-10	Passive	115	6F
MS1 - F0U - 06P6222	6	2-2-2	Passive	115	6F
MS1 - F0U - 06P6252	6	2-5-2	Passive	115	6F
MS1 - F0U - 06P6555	6	5-5-5	Passive	115	6F
MS1 - F0U - 06P6282	6	2-8-2	Passive	115	6F
MS1 - F0U - 06P62102	6	2-10-2	Passive	115	6F
MS1 - F0U - 06P6101010	6	10-10-10	Passive	115	6F
MS1 - F0U - 08P6222	8	2-2-2	Passive	115	6F
MS1 - F0U - 08P6252	8	2-5-2	Passive	115	6F
MS1 - F0U - 08P6555	8	5-5-5	Passive	115	6F
MS1 - F0U - 08P6282	8	2-8-2	Passive	115	6F
MS1 - F0U - 08P62102	8	2-10-2	Passive	115	6F
MS1 - F0U - 08P6101010	8	10-10-10	Passive	115	6F
MS1 - F0U - 10P6222	10	2-2-2	Passive	115	6F
MS1 - F0U - 10P6252	10	2-5-2	Passive	115	6F
MS1 - F0U - 10P6555	10	5-5-5	Passive	115	6F
MS1 - F0U - 10P6282	10	2-8-2	Passive	115	6F
MS1 - F0U - 10P62102	10	2-10-2	Passive	115	6F
MS1 - F0U - 10P6101010	10	10-10-10	Passive	115	6F

Steerable Electrophysiology Catheters With Reuse Functionality

G Curve Catheters



Compatible with
X-D10 Extension
Cable



Ordering Information


Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
MS1 - G0U - 04P6222	4	2-2-2	Passive	115	6F
MS1 - G0U - 04P6252	4	2-5-2	Passive	115	6F
MS1 - G0U - 04P6555	4	5-5-5	Passive	115	6F
MS1 - G0U - 04P6282	4	2-8-2	Passive	115	6F
MS1 - G0U - 04P62102	4	2-10-2	Passive	115	6F
MS1 - G0U - 04P6101010	4	10-10-10	Passive	115	6F
MS1 - G0U - 06P6222	6	2-2-2	Passive	115	6F
MS1 - G0U - 06P6252	6	2-5-2	Passive	115	6F
MS1 - G0U - 06P6555	6	5-5-5	Passive	115	6F
MS1 - G0U - 06P6282	6	2-8-2	Passive	115	6F
MS1 - G0U - 06P62102	6	2-10-2	Passive	115	6F
MS1 - G0U - 06P6101010	6	10-10-10	Passive	115	6F
MS1 - G0U - 08P6222	8	2-2-2	Passive	115	6F
MS1 - G0U - 08P6252	8	2-5-2	Passive	115	6F
MS1 - G0U - 08P6555	8	5-5-5	Passive	115	6F
MS1 - G0U - 08P6282	8	2-8-2	Passive	115	6F
MS1 - G0U - 08P62102	8	2-10-2	Passive	115	6F
MS1 - G0U - 08P6101010	8	10-10-10	Passive	115	6F
MS1 - G0U - 10P6222	10	2-2-2	Passive	115	6F
MS1 - G0U - 10P6252	10	2-5-2	Passive	115	6F
MS1 - G0U - 10P6555	10	5-5-5	Passive	115	6F
MS1 - G0U - 10P6282	10	2-8-2	Passive	115	6F
MS1 - G0U - 10P62102	10	2-10-2	Passive	115	6F
MS1 - G0U - 10P6101010	10	10-10-10	Passive	115	6F

Steerable Electrophysiology Catheters With Reuse Functionality

J Curve Catheters



Compatible with
X-D10 Extension
Cable



Ordering Information


Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
MS1 - J0U - 04P6222	4	2-2-2	Passive	115	6F
MS1 - J0U - 04P6252	4	2-5-2	Passive	115	6F
MS1 - J0U - 04P6555	4	5-5-5	Passive	115	6F
MS1 - J0U - 04P6282	4	2-8-2	Passive	115	6F
MS1 - J0U - 04P62102	4	2-10-2	Passive	115	6F
MS1 - J0U - 04P6101010	4	10-10-10	Passive	115	6F
MS1 - J0U - 06P6222	6	2-2-2	Passive	115	6F
MS1 - J0U - 06P6252	6	2-5-2	Passive	115	6F
MS1 - J0U - 06P6555	6	5-5-5	Passive	115	6F
MS1 - J0U - 06P6282	6	2-8-2	Passive	115	6F
MS1 - J0U - 06P62102	6	2-10-2	Passive	115	6F
MS1 - J0U - 06P6101010	6	10-10-10	Passive	115	6F
MS1 - J0U - 08P6222	8	2-2-2	Passive	115	6F
MS1 - J0U - 08P6252	8	2-5-2	Passive	115	6F
MS1 - J0U - 08P6555	8	5-5-5	Passive	115	6F
MS1 - J0U - 08P6282	8	2-8-2	Passive	115	6F
MS1 - J0U - 08P62102	8	2-10-2	Passive	115	6F
MS1 - J0U - 08P6101010	8	10-10-10	Passive	115	6F
MS1 - J0U - 10P6222	10	2-2-2	Passive	115	6F
MS1 - J0U - 10P6252	10	2-5-2	Passive	115	6F
MS1 - J0U - 10P6555	10	5-5-5	Passive	115	6F
MS1 - J0U - 10P6282	10	2-8-2	Passive	115	6F
MS1 - J0U - 10P62102	10	2-10-2	Passive	115	6F
MS1 - J0U - 10P6101010	10	10-10-10	Passive	115	6F

Steerable Electrophysiology Catheters With Reuse Functionality

P Curve Catheters



Compatible with
X-D10 Extension
Cable



Ordering Information


Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
MS1 - P0U - 04P6222	4	2-2-2	Passive	115	6F
MS1 - P0U - 04P6252	4	2-5-2	Passive	115	6F
MS1 - P0U - 04P6555	4	5-5-5	Passive	115	6F
MS1 - P0U - 04P6282	4	2-8-2	Passive	115	6F
MS1 - P0U - 04P62102	4	2-10-2	Passive	115	6F
MS1 - P0U - 04P6101010	4	10-10-10	Passive	115	6F
MS1 - P0U - 06P6222	6	2-2-2	Passive	115	6F
MS1 - P0U - 06P6252	6	2-5-2	Passive	115	6F
MS1 - P0U - 06P6555	6	5-5-5	Passive	115	6F
MS1 - P0U - 06P6282	6	2-8-2	Passive	115	6F
MS1 - P0U - 06P62102	6	2-10-2	Passive	115	6F
MS1 - P0U - 06P6101010	6	10-10-10	Passive	115	6F
MS1 - P0U - 08P6222	8	2-2-2	Passive	115	6F
MS1 - P0U - 08P6252	8	2-5-2	Passive	115	6F
MS1 - P0U - 08P6555	8	5-5-5	Passive	115	6F
MS1 - P0U - 08P6282	8	2-8-2	Passive	115	6F
MS1 - P0U - 08P62102	8	2-10-2	Passive	115	6F
MS1 - P0U - 08P6101010	8	10-10-10	Passive	115	6F
MS1 - P0U - 10P6222	10	2-2-2	Passive	115	6F
MS1 - P0U - 10P6252	10	2-5-2	Passive	115	6F
MS1 - P0U - 10P6555	10	5-5-5	Passive	115	6F
MS1 - P0U - 10P6282	10	2-8-2	Passive	115	6F
MS1 - P0U - 10P62102	10	2-10-2	Passive	115	6F
MS1 - P0U - 10P6101010	10	10-10-10	Passive	115	6F

Steerable Electrophysiology Catheters With Reuse Functionality

R Curve Catheters



Compatible with
X-D10 Extension
Cable



Ordering Information


Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
MS1 - R0U - 04P6222	4	2-2-2	Passive	115	6F
MS1 - R0U - 04P6252	4	2-5-2	Passive	115	6F
MS1 - R0U - 04P6555	4	5-5-5	Passive	115	6F
MS1 - R0U - 04P6282	4	2-8-2	Passive	115	6F
MS1 - R0U - 04P62102	4	2-10-2	Passive	115	6F
MS1 - R0U - 04P6101010	4	10-10-10	Passive	115	6F
MS1 - R0U - 06P6222	6	2-2-2	Passive	115	6F
MS1 - R0U - 06P6252	6	2-5-2	Passive	115	6F
MS1 - R0U - 06P6555	6	5-5-5	Passive	115	6F
MS1 - R0U - 06P6282	6	2-8-2	Passive	115	6F
MS1 - R0U - 06P62102	6	2-10-2	Passive	115	6F
MS1 - R0U - 06P6101010	6	10-10-10	Passive	115	6F
MS1 - R0U - 08P6222	8	2-2-2	Passive	115	6F
MS1 - R0U - 08P6252	8	2-5-2	Passive	115	6F
MS1 - R0U - 08P6555	8	5-5-5	Passive	115	6F
MS1 - R0U - 08P6282	8	2-8-2	Passive	115	6F
MS1 - R0U - 08P62102	8	2-10-2	Passive	115	6F
MS1 - R0U - 08P6101010	8	10-10-10	Passive	115	6F
MS1 - R0U - 10P6222	10	2-2-2	Passive	115	6F
MS1 - R0U - 10P6252	10	2-5-2	Passive	115	6F
MS1 - R0U - 10P6555	10	5-5-5	Passive	115	6F
MS1 - R0U - 10P6282	10	2-8-2	Passive	115	6F
MS1 - R0U - 10P62102	10	2-10-2	Passive	115	6F
MS1 - R0U - 10P6101010	10	10-10-10	Passive	115	6F

Steerable Electrophysiology Catheters With Reuse Functionality

P Curve Catheters



Compatible with
X-D20 Extension
Cable



Ordering Information


Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
MS2 - P0U - 20A6242	20	2-4-2	Active	115	6F
MS2 - P0U - 20A6282	20	2-8-2	Active	115	6F
MS2 - P0U - 20A62102	20	2-10-2	Active	115	6F
MS2 - P0U - 20A62122	20	2-12-2	Active	115	6F
MS2 - P0U - 20A62132	20	2-13-2	Active	115	6F
MS2 - P0U - 20P6242	20	2-4-2	Passive	115	6F
MS2 - P0U - 20P6282	20	2-8-2	Passive	115	6F
MS2 - P0U - 20P62102	20	2-10-2	Passive	115	6F
MS2 - P0U - 20P62122	20	2-12-2	Passive	115	6F
MS2 - P0U - 20P62132	20	2-13-2	Passive	115	6F
MS2 - P0U - 20P620282	20	20-2-8-2	Passive	115	6F

Steerable Electrophysiology Catheters With Reuse Functionality

Q Curve Catheters



Compatible with
X-D20 Extension
Cable



Ordering Information


Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
MS2 - Q0U - 20A6242	20	2-4-2	Active	115	6F
MS2 - Q0U - 20A6282	20	2-8-2	Active	115	6F
MS2 - Q0U - 20A62102	20	2-10-2	Active	115	6F
MS2 - Q0U - 20A62122	20	2-12-2	Active	115	6F
MS2 - Q0U - 20A62132	20	2-13-2	Active	115	6F
MS2 - Q0U - 20P6242	20	2-4-2	Passive	115	6F
MS2 - Q0U - 20P6282	20	2-8-2	Passive	115	6F
MS2 - Q0U - 20P62102	20	2-10-2	Passive	115	6F
MS2 - Q0U - 20P62122	20	2-12-2	Passive	115	6F
MS2 - Q0U - 20P62132	20	2-13-2	Passive	115	6F
MS2 - Q0U - 20P620282	20	20-2-8-2	Passive	115	6F

Steerable Electrophysiology Catheters With Reuse Functionality

R Curve Catheters



Compatible with
X-D20 Extension
Cable




Ordering Information

Product Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
MS2 - R0U - 20A6242	20	2-4-2	Active	115	6F
MS2 - R0U - 20A6282	20	2-8-2	Active	115	6F
MS2 - R0U - 20A62102	20	2-10-2	Active	115	6F
MS2 - R0U - 20A62122	20	2-12-2	Active	115	6F
MS2 - R0U - 20A62132	20	2-13-2	Active	115	6F
MS2 - R0U - 20P6242	20	2-4-2	Passive	115	6F
MS2 - R0U - 20P6282	20	2-8-2	Passive	115	6F
MS2 - R0U - 20P62102	20	2-10-2	Passive	115	6F
MS2 - R0U - 20P62122	20	2-12-2	Passive	115	6F
MS2 - R0U - 20P62132	20	2-13-2	Passive	115	6F
MS2 - R0U - 20P620282	20	20-2-8-2	Passive	115	6F

Disposable Sheaths for Steerable Catheters With Reuse Functionality

Sheaths are compatible with MS1 Catheters

MS2
Catheter
Compatible




Ordering Information

Sheath Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
04P6222	4	2-2-2	Passive	115	6F
04P6252	4	2-5-2	Passive	115	6F
04P6555	4	5-5-5	Passive	115	6F
04P6282	4	2-8-2	Passive	115	6F
04P62102	4	2-10-2	Passive	115	6F
04P6101010	4	10-10-10	Passive	115	6F
06P6222	6	2-2-2	Passive	115	6F
06P6252	6	2-5-2	Passive	115	6F
06P6555	6	5-5-5	Passive	115	6F
06P6282	6	2-8-2	Passive	115	6F
06P62102	6	2-10-2	Passive	115	6F
06P6101010	6	10-10-10	Passive	115	6F
08P6222	8	2-2-2	Passive	115	6F
08P6252	8	2-5-2	Passive	115	6F
08P6555	8	5-5-5	Passive	115	6F
08P6282	8	2-8-2	Passive	115	6F
08P62102	8	2-10-2	Passive	115	6F
08P6101010	8	10-10-10	Passive	115	6F
10P6222	10	2-2-2	Passive	115	6F
10P6252	10	2-5-2	Passive	115	6F
10P6555	10	5-5-5	Passive	115	6F
10P6282	10	2-8-2	Passive	115	6F
10P62102	10	2-10-2	Passive	115	6F
10P6101010	10	10-10-10	Passive	115	6F

Disposable Sheaths for Steerable Catheters With Reuse Functionality

Sheaths are compatible with MS2 Catheters




Ordering Information

Sheath Reference Code	No. of Electrodes	Electrode Spacing (mm)	Distal Type	Working Length (cm)	French Size
20A6242	20	2-4-2	Active	115	6F
20A6282	20	2-8-2	Active	115	6F
20A62102	20	2-10-2	Active	115	6F
20A62122	20	2-12-2	Active	115	6F
20A62132	20	2-13-2	Active	115	6F
20P6242	20	2-4-2	Passive	115	6F
20P6282	20	2-8-2	Passive	115	6F
20P62102	20	2-10-2	Passive	115	6F
20P62122	20	2-12-2	Passive	115	6F
20P62132	20	2-13-2	Passive	115	6F
20P620282	20	20-2-8-2	Passive	115	6F

Loop Steerable Electrophysiology Catheter

STEP 1

In the event of a difficult catheter retraction physicians may release the loop by separating the electrode sheath from the catheter handle.




STEP 2

The loop loses its support and unravels, becoming very flexible. The catheter can now be withdrawn.

Loop Steerable Electrophysiology Catheter

All loop catheters have a passive tip



Ordering Information


Reference Code	Curve Type	Loop Size (mm)	Electrodes	Electrode Spacing	Shaft Size
SC2 - B2515U - 10P7888	B	25-15	10	8-8-8	7F
SC1 - B15U - 10P7555	B	15	10	5-5-5	7F
SC1 - B20U - 10P7666	B	20	10	6-6-6	7F
SC1 - B25U - 10P7888	B	25	10	8-8-8	7F
SC1 - D15U - 10P7555	D	15	10	5-5-5	7F
SC1 - D20U - 10P7666	D	20	10	6-6-6	7F
SC1 - D25U - 10P7888	D	25	10	8-8-8	7F

Loop Steerable Electrophysiology Catheter

All loop catheters have a passive tip



Compatible with
X-D20 Extension
Cable




Ordering Information


Reference Code	Curve Type	Loop Size (mm)	Electrodes	Electrode Spacing	Shaft Size
SC2 - B2515U - 20P7262	B	25-15	20	2-6-2	7F
SC2 - B15U - 14P7333	B	15	14	3-3-3	7F
SC2 - B20U - 14P7444	B	20	14	4-4-4	7F
SC2 - B30U - 14P7777	B	30	14	7-7-7	7F
SC2 - B15U - 20P7232	B	15	20	2-3-2	7F
SC2 - B20U - 20P7242	B	20	20	2-4-2	7F
SC2 - B25U - 20P7262	B	25	20	2-6-2	7F
SC2 - B30U - 20P7272	B	30	20	2-7-2	7F
SC2 - D15U - 14P7333	D	15	14	3-3-3	7F
SC2 - D20U - 14P7444	D	20	14	4-4-4	7F
SC2 - D30U - 14P7777	D	30	14	7-7-7	7F
SC2 - D15U - 20P7232	D	15	20	2-3-2	7F
SC2 - D20U - 20P7242	D	20	20	2-4-2	7F
SC2 - D25U - 20P7262	D	25	20	2-6-2	7F
SC2 - D30U - 20P7272	D	30	20	2-7-2	7F

Extension Cables

10 Pin Extension Cable



20 Pin Extension Cable



Reference Code	No. of Connections	Description
X-D10	10	Extension cable for up to 10 electrode diagnostic catheters
X-D20	20	Extension cable for 11 - 20 electrode diagnostic catheters

Sterilisation Trays

Reference Code	Accessory Type	Description
A-T01	Tray	Sterilisation Tray
A-T02	Tray	Washing and Sterilisation Tray

